

**Shiga Toxin-producing *E. coli* (STEC),
O157 and Non-O157: Efforts in
Wisconsin to
Improve Detection**

Dave Warshauer, PhD and Tim Monson, MS
Communicable Disease Division
Wisconsin State Laboratory of Hygiene
(608)265-9115/ (608)263-3421



WSLH Audio Conference Program
March 17, 2004



Objectives

- Describe the epidemiology, pathogenesis and clinical manifestations of STEC disease
- Describe various testing protocols for the detection of STEC infection
- Describe testing available at WSLH for the detection of STEC infection
- Describe a state-wide collaborative testing program for STEC detection and surveillance

WSLH

All *E. coli* are not created equal

- Several *E. coli* groups can cause outbreaks of severe diarrhea:

STEC, EHEC	enterohemorrhagic (ie: O157:H7)
ETEC	enterotoxigenic
EPEC	enteropathogenic
EAEC	enteroadherent
EaggEC	enteroaggregative
EIEC	enteroinvasive
DAEC	diffuse adherant

WSLH

History

- 1885----E. coli first described
- 1982----E. coli O157:H7 first recognized as a human pathogen
- 1985----Karmali linked toxin-producing E.coli to hemolytic uremic syndrome (HUS)
- 1993----Outbreak in four western states affecting 500+ persons, 56 HUS, 4 deaths
 - Led to recommendations for screening

WSLH

STEC

- A shiga toxin-producing E. coli (STEC)
 - AKA verotoxin, verocytotoxin, shiga-like toxin
 - Stx1 and Stx2
 - Similar to shiga toxin produced by *S. dysenteriae*
- Possess other virulence factors
 - Attaching and effacing gene (eae)
 - Enterohemolysin gene (E-hly)
- Other E. coli serotypes can produce shiga toxin
 - O26:H11 O145:NM O45:H2
 - O111:H8 O118:H2 O104:H21
 - O113:H21 O121:H19
 - >100 serotypes

WSLH

E. coli O157:H7

- O157:H7 and O157:NM
 - 73,000 cases of illness each year in U.S.
 - 60 deaths
 - 80% of cases of hemolytic uremic syndrome (HUS)
- Prevalence in Wisconsin----CHOW
 - Kehl et. al (JCM 35:2051-54, 1997)
 - E. coli O157:H7-----1.2%
 - Non-O157 STEC-----0.3%
 - Salmonella-----2.8%
 - Shigella-----3.1%
 - Campylobacter-----0.9%

WSLH

Cases of enteric pathogens reported in Wisconsin 2002

- Campylobacter-----1159
- Salmonella----- 876
- E. coli O157:H7----- 293
- Shigella----- 178
- non-O157:H7 STEC----- 5

WSLH

Surveillance

- 1993--CSTE recommended that infections be made reportable by all states
- 2000--Made a reportable disease in Wisconsin
- Nationally reportable disease
- Part of the CDC PulseNet Program
 - Established 1996

WSLH

Clinical Aspects

- Incubation period
 - Typically 2-4 days, range 1-7
 - HUS in 5-10% of infected
 - presents 2-10 days following the initial infection
- Initial symptoms
 - Diarrhea--mild, watery, bloody
 - Fecal leukocytes in >70%
 - Abdominal cramps
 - Vomiting in about 50% of patients

WSLH

Hemolytic Uremic Syndrome

- More frequent in young (<5 years) and the elderly (>75 years)
 - Unclear why, but may be related to lack of immunity
- Triad of characteristics
 - renal failure
 - thrombocytopenia
 - hemolytic anemia

WSLH

Outcome of HUS

- 5-10% mortality
 - Multi-organ failure---CNS, pancreas, heart, and others
- Long-term complication rate of 30-50%
 - Chronic renal failure
 - proteinuria
 - Hypertension
- Stroke----3-4%

WSLH

Other complications

- Thrombotic thrombocytopenic purpura (TTP)
much like HUS except:
 - occurs in a more adult population
 - associated with neurological complications
- Thrombotic strokes
- Blindness
- Coma
- Cardiac failure
- Gut perforation, gangrene

WSLH

Pathogenesis

- Colonize portions of the lower intestine
 - Capacity to adhere, attach, and efface epithelial cells important
- Production of Shiga-toxin (STEC)
 - Cross intestinal epithelium, reaching sub-epithelial tissues
- Endothelial cells lining small blood vessels are final targets for the toxins
 - kidney, brain, and other organs

WSLH

Pathogenesis (2)

- Formation of aggregates of platelets within blood vessels with generation of thrombotic microangiopathy

WSLH

Transmission

- Cattle primary reservoir
 - Prevalence can reach 63-100% in feedlots
 - Part of normal intestinal flora
 - No clinical illness
- Contamination of meat during slaughter
 - Spread through grinding and mixing
- Transmission to humans
 - Consumption of undercooked beef or beef products
 - Cross contamination of other foods

WSLH

Transmission (2)

- Cases or outbreaks associated with
 - Beef, salami, sausage, other meat products
 - unpasteurized milk
 - unpasteurized apple juice/cider
 - lettuce
 - alfalfa and radish sprouts
 - water---unchlorinated drinking and recreational
- Person-to-person transmission
 - Shedding—most <1 week; some >3 weeks

WSLH

Outbreaks

- 52% Cattle-related food
 - 45% ground beef
- 14% Fruits and vegetables
- 12% Water
- 2% Milk
- 16% Direct transmission
 - person-to-person
 - cattle-to-person

WSLH

Important Outbreaks of non-O157:H7 STEC

- Australia, 1995 (see Paton *et al.*, JCM 34:1622-7)
 - >160 cases
 - 22 HUS, 4 TTP
 - O111:H-, stx1+, 2+, eae+ (PCR)
 - Mettwurst (dry, fermented sausage) was infection source
- Texas, 1999 (see MMWR 49:321-4)
 - 58 cases,
 - 2 HUS
 - O111:H8, stx1+, 2+ (PCR and EIA)
 - Ice in drinking water implicated as source
- Outbreaks of sorbitol-fermenting O157:H- in Europe (see Ammon *et al.*, JID 179: 1274-7)

WSLH

Wisconsin Outbreaks
E. coli O157:H7

- Ozaukee County Fair, Cedarburg, August 2001
 - E. coli O157:H7
 - Associated with contact with cattle
 - 25 lab confirmed cases
 - 34 probable
 - No HUS or fatalities
 - Led to state guidelines for exhibit of animals

WSLH

Wisconsin Outbreaks (2)
E. coli O157:H7

- UW-Madison Pancake Breakfast, October 2001
 - Held in the stock pavilion on campus
 - 34 ill, 16 confirmed by culture
 - 1 child with HUS
 - No food implicated
 - Pavilion houses animals
 - 3 days previous, cows housed in a pen in the pavilion where the foodline was located
 - Outbreak strain isolated from environment and one of the cows

WSLH

Wisconsin Outbreaks (3)
E. coli O157:H7

- 2 Multi-county, multi-state outbreaks, Summer 2001
 - 9 states involved in one outbreak, 13 in other
 - 2 different strains of E. coli O157:H7
 - 54 lab confirmed WI cases
 - 44 cases from other states
 - 1 HUS, no fatalities
 - Associated with consumption of ground beef
 - One outbreak strain isolated from ground beef
 - Led to meat recall

WSLH

Proportion of STEC isolates that were O157 vs non-O157, diarrheal stools

State	STEC	O157	non-O157
MN	134	74 (55%)	60 (45%)
CT	32	19 (59%)	13 (41%)
NM	35	14 (40%)	21 (60%)
Neb	14	7 (50%)	7 (50%)
Total	215	114 (53%)	101 (47%)

- MN – Besser et al., unpublished
- CT - Fiorentino et al., ASM 2001 (Abstract C166)
- NM – Nims et al. ASM 2001 (Abstract C165)
- Neb – Fey et al. EID 6:530-533, 2000

WSLH

CDC Recommendations
(MMWR, 4/21/00)

- Clinicians should inform health departments about clusters of bloody diarrhea and HUS
- Clinical Labs should screen bloody stools for STEC, and attempt to isolate STEC from positive stools
- Positive isolates should be referred to public health laboratories for serotyping and further characterization.
- States should consider adding STEC infection to their notifiable disease lists.

WSLH



Tim Monson


WSLH

Laboratory Diagnosis of STEC: Culture for *E. coli* O157:H7

- Serotype O157-Specific Culture Methods
 - Selective and differential media
 - **Sorbitol MacConkey** agar plates (SMAC)
 - **Cefixime-Tellurite Sorbitol MacConkey** Agar (CT-SMAC)
- Chromagars:
 - Rainbow Agar O157
 - ChromAgar O157


WSLH

Chromagar O157



WSLH

Rainbow Agar



WSLH

**Laboratory Diagnosis of STEC:
Culture for *E. coli* O157:H7**

- *E. coli* O157:H7 either fail to ferment sorbitol or exhibit delayed sorbitol fermentation (rare exceptions occur)
- Non-O157 bacteria that grow on SMAC may also be sorbitol negative
- Must confirm using specific agglutination assays for O157 and H7 antigens

WSLH

**Laboratory Diagnosis of STEC:
Culture for *E. coli* O157:H7**

- Biochemical confirmation of *E. coli* O157:H7
 - **Methylumbelliferyl- β -D-glucuronide (MUG)**
Test to test for the presence of β -glucuronidase enzyme (O157:H7 are typically MUG negative)
 - Common biochemical profile for O157:H7 strains

WSLH

**Laboratory Diagnosis of STEC:
Disadvantages of Conventional
Testing for *E. coli* O157:H7**

- Sorbitol-fermenting *E. coli* O157:H7 strains
- Low numbers of O157:H7 may not be isolated
- Sensitivity of culturing for O157:H7 strains, using SMAC plate only, is estimated to be only 50-60%
- Non-O157:H7 STEC strains will not be detected


WSLH

Laboratory Diagnosis of STEC: Enhanced Culture Method

- Immunomagnetic Separation (IMS)
 - Metallic beads coated with anti-O (LPS) antibody (anti-O157, anti-O26, anti-O111, etc)
 - Mix beads with stool in enrichment broth
 - Beads specifically bind bacterial LPS
 - Immobilize beads with a magnetic field
 - Wash then plant beads on plate media
 - May increase sensitivity of culture 100-fold

WSLH

Immunomagnetic Separation -Scanning electron-micrograph



WSLH

Laboratory Diagnosis of STEC: Non-culture Methods

- Toxin Assays
 - Vero cells commonly used
- EIA
 - Meridian Premier EHEC
 - Alexon Trend ProSpecT Shiga Toxin E. coli (STEC)
- RPLA
 - Oxoid VTEC-RPLA

WSLH


ImmunoCard STAT!
E. coli O157:H7



WSLH

Reversed Passive Latex
Agglutination Assays

- Latex particles sensitized with anti-serum reactive with Shiga toxins (verocytotoxins) 1 and 2
- Agglutination of latex particles and toxin results in a lattice formation at the bottom of the V-shaped microtiter plate well



WSLH

Laboratory Diagnosis of STEC:
Non-serotype specific testing

- Non-serotype specific assays involve:
 - (1) Detection of Shiga toxins (Toxin assays)
 - (2) Detection of specific virulence genes associated with STEC disease (Genetic assays)
- Non-serotype specific assays will detect non-O157:H7 STEC strains as well as O157:H7 strains

WSLH

Laboratory Diagnosis of STEC:
Toxin Assays

- Cytotoxicity in cell culture
 - Vero cells commonly used
- Immunoassays
 - Commercially available--- Alexon-Trend, Meridian Biosciences, Oxoid
 - Test stools directly--- (70-80% sensitive)
 - Incubate in enrichment broth (GN or MAC) --- (approach 100% sensitivity)

WSLH

Isolation and Confirmation of STEC:
WSLH Testing Algorithm


- Stool specimen submitted
 - (1) Plate to MAC, SMAC and CT-SMAC to rule out *E. coli* O157:H7
 - (2) GN enrichment broth set up for Shiga toxin EIA testing
 - (3) Shiga toxin EIA assay performed if culture negative for *E. coli* O157:H7
 - (4) Enrichment broths positive for Shiga toxin are plated for PCR testing and isolation of STEC

WSLH


Isolation and Confirmation of STEC: WSLH Testing Algorithm

- GN Enrichment broth (Shiga toxin positive) or isolate submitted
 - (1) Plated to MAC, SMAC and CT-SMAC to rule out *E. coli* O157:H7
 - (2) Sweep of growth from MAC or SMAC plate tested by PCR for the presence of Shiga toxin genes (*stx1* and *stx2*)
 - (3) If positive by PCR, individual colonies tested to find STEC


WSLH



- ### E. coli O157 and STEC Virulence Genes
- *stx1* Shiga toxin 1 (SLT-1/ VT 1)
 - *stx2* Shiga toxin 2 (SLT-2/ VT 2)
 - *eaeA* *E. coli* attaching and effacing gene (Intimin)
 - *uidA* O157 allele of beta glucuronidase
 - *E-hly* Enterohemolysin (found on plasmid present in *E. coli* O157:H7 and other EHEC)
- WSLH



Real-Time PCR Assay



WSLH

Laboratory Diagnosis of STEC: Serotyping


- Slide Agglutination at WSLH
 - EPEC/VTEC Pool USA
 - O45 and O121
 - EPEC/VTEC Pool 1
 - O26, O103, O111, O145 and O157
- STEC sent to CDC for H antigen serotyping

WSLH

STEC Serotypes isolated in WI

- O157:H7 (many)
- O28ac:H25
- O45:H(pending)
- O45:H2
- O103:H2 (4)
- O111:NM
- O157:NM (many)
- O145:NM (2)
- O?:H28
- O?:NM
- unnamed isolate
- O and H pending

WSLH



Wisconsin SHIELD Program

- **Shiga Toxin Epidemiology and Laboratory Detection Program**
- Surveillance program to determine the prevalence of STEC disease in WI and identify/serotype the agents involved
- Participating laboratories are sending stool specimens, bacterial isolates and enrichment broths from patients that meet the program criteria

WSLH

Wisconsin SHIELD Program: Submission Criteria

- Specimens from patients with a history of diarrheal illness that also:
 - (1) Are suspected or confirmed HUS cases
 - (2) Are suspected or confirmed of containing blood
 - (3) Are from pediatric patients
 - (4) Are known or suspected contacts of confirmed STEC cases or outbreaks

WSLH

WSLH WEPS Program



WSLH

WEPS Program- STEC

- **Wisconsin Enteric Pathogens Surveillance Program**- requests submission of STEC specimens and other enteric pathogens to WSLH for PFGE analysis and antimicrobial resistance monitoring
- Shiga toxin EIA positive enrichment broths or resulting STEC isolate
- Confirmed or suspect STEC Isolates requested
 - *E. coli* O157:H7
 - Non-O157 STECS

WSLH

Submission of Specimens

- Dunham Express courier service is available for delivery of SHIELD or WEPS Program specimens to WSLH
 - Account number 7271
 - Call 1-800-236-7127

WSLH

References

- Kehl, Sue C. 2002. Role of the laboratory in the diagnosis of enterohemorrhagic *Escherichia coli* infections. *J. Clin. Microbiol.* 40(8):2711-2715.
- Paton, James C and Adrienne W Paton. 1998. Pathogenesis and diagnosis of Shiga toxin-producing *Escherichia coli* infections. *Clin. Micro. Rev.* 11(3):450-479.

_____ WSLH

References- continued

- Park, C.H., H.J. Kim and D.L. Hixon. 2002. Importance of testing stool specimens for shiga toxins. *J. Clin. Microbiol.* 40(9):3542-3543.
- Strockbine, N.A., et al. Overview of detection and subtyping methods. *In: J.B. Kapler and A.D. O'Brien ed. Escherichia coli O157:H7 and other Shiga toxin-producing E. coli strains.* Washington, D.C.: ASM Press; 1998: 331-356

_____ WSLH

Contact Information

- Dave Warshauer, PhD. ABMM
 - Assistant Director- Comm. Disease Division
 - (608)265-9115
 - warshadmin@mail.slh.wisc.edu
- Tim Monson, MS
 - Advanced Microbiologist- Food- and Waterborne Disease Program Coordinator
 - (608)263-3421
 - mony@mail.slh.wisc.edu

_____ WSLH
