

2024 Updates to CLSI M100



Erik Munson
Marquette University
Wisconsin Clinical Laboratory Network
Laboratory Technical Advisory Group

The presenter states no conflict of interest and has no financial relationship to disclose relevant to the content of this presentation.

OUTLINE

I. Quick discussion(s) relative to major revisions

II. Objectives of webinar


Describe significant changes relevant to pre-existing antimicrobial susceptibility breakpoints...

Describe significant changes relevant to antimicrobial susceptibility testing methodology...

Identify (new) organism/antimicrobial combinations for which susceptibility breakpoints now exist...

as outlined in the CLSI M100-Ed34 document.

clsi.org/m100, then scroll down quite a bit



“CLSI FREE”...then a couple of clicks

Free Resources From CLSI

https://clsi.org/all-free-resources/

New editions of CLSI M100, M02, and M07 now available with discount package pricing.


eCLIPSE eLearning Exchange Support Meetings Sign in

CLSI

Shop Membership Participate Standards Global Training About

Access Our Free Resources

Free Microbiology Resources



Micro Free

CLSI Micro Free
M100, M23, M23 Supplements, M27M44S, and M45

Micro Free (formerly M100 Free) includes read-only access to critical microbiology standards. These standards are essential for better patient outcomes and reducing antimicrobial resistance. CLSI is committed to making resources broadly available to users around the globe to save lives and improve patient health.

[Access Now →](#)

CLICK HERE



VET015
Performance Standards for Antimicrobial Disk and Dilution Susceptibility Tests for Bacteria Isolated From Animals

Quickly reference the most trusted AST veterinary breakpoint tables as a convenient, complimentary supplement to the AST [VET01 document](#).

[Access Now →](#)



CLSI AST Rationale Documents
Package of Rationale Documents

Providing the scientific reasons behind breakpoint decisions



CLSI Using M100-Ed32
Using M100-Ed32

BIG ADVERTISEMENT...

CLSI eClipse Ultimate Access - Po x

https://em100.edaptivedocs.net/Login.aspx?_ga=2.174075298.1086753664.1709112633-1991978545.1705003341

Email:

Login

CLSI MICRO FREE

CLSI MicroFree enables free access to the most trusted and updated information on Antimicrobial Susceptibility Testing.

This platform is offered as a public health service. For more complete standards documents, training, and related products (including digital, print, and subscription versions of these and other AST documents, visit clsi.org/standards).

© 2024 CLSI | © 2024 Edaptive Technologies LLC

Powered By: Edaptive Platform

NEW

AST FREE PORTAL

CLSI

Sign-Out

Welcome to the AST Free Portal!

Enter Keywords Search

My Library

- [CLSI M100 ED34:2024](#)
- [CLSI M45 ED3:2016](#)
- [CLSI M27M44S ED3:2022](#)
- [CLSI M23S3 ED1:2023](#)
- [CLSI M23S2 ED1:2021](#)
- [CLSI M23S ED1:2020](#)
- [CLSI M23 ED6:2023](#)

Bulletin Board

Document Updates

2024-02-28	CLSI M100 ED34:2024
2023-07-28	CLSI M23 ED6:2023
2023-06-15	CLSI M23S3 ED1:2023

[See more ...](#)

Message

2015-11-10 **Welcome to the CLSI MicroFree Portal**


This platform provides free access the most trusted and updated information on Antimicrobial

M27M44S involves broth microdilution, disk diffusion susceptibility testing of yeast

TABLE OF CONTENTS

The screenshot shows a web browser window with the following elements:

- Browser Tab:** EM100 Connect - CLSI M100 ED34
- Address Bar:** <https://em100.edaptivedocs.net/GetDoc.aspx?doc=CLSI%20M100%20ED34:2024&sbssok=CLSI%20M100%20ED34:2024%20SECTIO...>
- Header:** CLSI logo on the left, "Sign-Out" and "Home" links on the right.
- Left Navigation Panel:**
 - CLSI M100 ED34:2024
 - Go To: Top | Bottom
 - Table of Contents (highlighted with a blue bar and a minus sign)
 - Table 2A-1. Zone Diameter and MIC Breakpoints for Enterobacterales (excluding *Salmonella/Shigella*)
 - Table 2A-2. Zone Diameter and MIC Breakpoints for *Salmonella* and *Shigella* spp.
 - Table 2B-1. Zone Diameter and MIC Breakpoints for *Pseudomonas aeruginosa*
 - Table 2B-2. Zone Diameter and MIC Breakpoints for *Acinetobacter* spp.
 - Table 2B-3. MIC Breakpoints for *Burkholderia cepacia* complex
 - Table 2B-4. Zone Diameter and MIC Breakpoints for *Stenotrophomonas maltophilia*
 - Table 2B-5. MIC
 - Related Information (highlighted with a blue bar and a plus sign)
- Main Content Area:**
 - CLSI M100-ED34:2024 Performance Standards for Antimicrobial Susceptibility Testing, 34th Edition
 - Search within this Document (input field)
 - < Previous | Next >
 - Introduction**
 - CLSI M100 includes updated tables for the Clinical and Laboratory Standards Institute antimicrobial susceptibility testing standards CLSI M02, M07, and M11.
 - A CLSI supplement for global application.
 - CLSI M100-Ed34**
 - February 2024**
 - Replaces CLSI M100-Ed33**
 - James S. Lewis II, PharmD, FIDSA
 - Amy J. Mathers, MD, D(ABMM)
 - April M. Bobenchik, PhD, D(ABMM)
 - Alexandra Lynn Bryson, PhD, D(ABMM)
 - Shelley Campeau, PhD, D(ABMM)
 - Sharon K. Cullen, BS, RAC
 - Tanis Dingle, PhD, D(ABMM), FCCM
 - Marcelo F. Galas, BSc
 - Romney M. Humphries, PhD, D(ABMM), FIDSA
 - Thomas J. Kirn, Jr., MD, PhD



Pseudomonas aeruginosa

EM100 Connect - CLSI M100 ED34:2024

https://em100.edaptivedocs.net/GetDoc.aspx?doc=CLSI%20M100%20ED34:2024&sbssok=CLSI%20M100%20ED34:2024%20TABLE%...

CLSI

Sign-Out

Home

CLSI M100 ED34:2024

Go To: Top | Bottom

Table of Contents

- Introduction
- Abstract
- Committee Membership
- Acknowledgment
- Overview of Changes
- CLSI Breakpoint Additions Since 2010
- CLSI Breakpoint Revisions Since 2010
- CLSI Archived Resources
- Summary of CLSI Processes for Establishing Breakpoints and QC Ranges
- CLSI Reference Methods vs Commercial Methods and CLSI vs US Food and Drug Administration Breakpoints
- CLSI Subcommittee on Antimicrobial

Related Information

Table 2B-1. Zone Diameter and MIC Breakpoints for *Pseudomonas aeruginosa*

Testing Conditions	Routine QC Recommendations (see Tables 4A-1 and 5A-1 for acceptable QC ranges)
Medium: Disk diffusion: MHA Broth dilution: CAMHB; iron-depleted CAMHB for cefiderocol (see Appendix H) ¹ Agar dilution: MHA	<i>P. aeruginosa</i> ATCC ^{®a} 27853
Inoculum: Broth culture method or colony suspension, equivalent to a 0.5 McFarland standard; positive blood culture broth for select antimicrobial agents with disk diffusion (see general comment [6])	Refer to Tables 4A-2 and 5A-2 to select strains for routine QC of β -lactam combination agents.
Incubation: 35°C \pm 2°C; ambient air Disk diffusion: 16–18 hours Dilution methods: 16–20 hours	When a commercial test system is used for susceptibility testing, refer to the manufacturer's instructions for QC test recommendations and QC ranges.

Refer to Tables 3B, 3C, 3E, 3F-1, and 3F-3 for additional testing recommendations, reporting suggestions, and QC.

General Comments

- (1) Refer to Table 1B-1 for antimicrobial agents that should be considered for testing and reporting by microbiology laboratories.
- (2) For disk diffusion, test a maximum of 12 disks on a 150-mm plate and no more than 6 disks on a 100-mm plate; disks should be placed no less than 24 mm apart, center to center (see CLSI M02²). Each zone diameter should be clearly

Pseudomonas aeruginosa



Sign-Out

Home

CLSI M100 ED34:2024

Go To: Top | Bottom

Table of Contents

- Introduction
- Abstract
- Committee Membership
- Acknowledgment
- Overview of Changes
- CLSI Breakpoint Additions Since 2010
- CLSI Breakpoint Revisions Since 2010
- CLSI Archived Resources
- Summary of CLSI Processes for Establishing Breakpoints and QC Ranges
- CLSI Reference Methods vs Commercial Methods and CLSI vs US Food and Drug Administration Breakpoints
- CLSI Subcommittee on Antimicrobial

Related Information

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, µg/mL			Comments
		S	I	R	S	I	R	
PENICILLINS								
Piperacillin*	100 µg	≥ 22	18-21 [^]	≤ 17	≤ 16	32 [^]	≥ 64	
β-LACTAM COMBINATION AGENTS								
(7) Organisms that test susceptible to the β-lactam agent alone are also considered susceptible to the β-lactam combination agent. However, organisms that test susceptible to the β-lactam combination agent cannot be assumed to be susceptible to the β-lactam agent alone. Similarly, organisms that test intermediate or resistant to the β-lactam agent alone may be susceptible to the β-lactam combination agent.								
Piperacillin-tazobactam	100/10 µg	≥ 22	18-21	≤ 17	≤ 16/4	32/4	≥ 64/4	(8) Breakpoints for intermediate are only to provide a buffer zone to prevent small uncontrolled technical factors from causing major discrepancies in interpretation.
Ceftazidime-avibactam	30/20 µg	≥ 21	-	≤ 20	≤ 8/4	-	≥ 16/4	
Ceftolozane-tazobactam	30/10 µg	≥ 21	17-20 [^]	≤ 16	≤ 4/4	8/4 [^]	≥ 16/4	
Imipenem-relebactam	10/25 µg	≥ 23	20-22 [^]	≤ 19	≤ 2/4	4/4 [^]	≥ 8/4	
Ticarcillin-clavulanate*	75/10 µg	≥ 24	16-23 [^]	≤ 15	≤ 16/2	32/2-64/2 [^]	≥ 128/2	
CEPHEMS (PARENTERAL) (Including cephalosporins I, II, III, and IV. Please refer to Glossary I.)								
Ceftazidime	30 µg	≥ 18	15-17 [^]	≤ 14	≤ 8	16 [^]	≥ 32	
Cefepime	30 µg	≥ 18	15-17 [^]	≤ 14	≤ 8	16 [^]	≥ 32	

RELATED INFORMATION

The screenshot shows a web browser displaying the CLSI M100 ED34:2024 Performance Standards for Antimicrobial Susceptibility Testing, 34th Edition. The page features a blue header with the CLSI logo and navigation links for 'Sign-Out' and 'Home'. A sidebar on the left contains a 'Table of Contents' and 'Related Information' section. The 'Related Information' section lists various CLSI standards and documents, with 'CLSI M45' highlighted by a black arrow. The main content area includes a search bar, navigation arrows, and a list of related information items.

CLSI M100 ED34:2024
Go To: [Top](#) | [Bottom](#)
[Table of Contents](#) +
[Related Information](#) -

Document Status:
Active

References:

- CLSI M02
- CLSI M07
- CLSI M11
- [CLSI M23](#)
- CLSI M39
- [CLSI M45](#)
- CLSI M52
- IDSA 2023
- ISO 20776-1

Referenced by:

- [CLSI M23S ED1:2020](#)

CLSI M100-ED34:2024 Performance Standards for Antimicrobial Susceptibility Testing, 34th Edition

Search within this Document

< Previous | Next >

Introduction

CLSI M100 includes updated tables for the Clinical and Laboratory Standards Institute antimicrobial susceptibility testing standards CLSI M02, M07, and M11.


A CLSI supplement for global application.

CLSI M100-Ed34
February 2024
Replaces CLSI M100-Ed33

James S. Lewis II, PharmD, FIDSA
Amy J. Mathers, MD, D(ABMM)
April M. Bobenchik, PhD, D(ABMM)
Alexandra Lynn Bryson, PhD, D(ABMM)
Shelley Campeau, PhD, D(ABMM)
Sharon K. Cullen, BS, RAC
Tanis Dingle, PhD, D(ABMM), FCCM
Marcelo F. Galas, BSc
Romney M. Humphries, PhD, D(ABMM), FIDSA
Thomas J. Kirn, Jr., MD, PhD

M02
M07
M11
M23
M39
M45
M52

Disk diffusion
Broth microdilution
Anaerobes
Quality control
Antibiogram
Infrequent bacteria
Commercial verification



FOR EXAMPLE, DIRECT LINK TO M45

The screenshot shows a web browser window with the URL <https://em100.edaptivedocs.net/GetDoc.aspx?doc=CLSI%20M45%20ED3:2016&format=HTM>. The page features the CLSI logo in the top left and navigation links for 'Sign-Out' and 'Home' in the top right. A left sidebar contains a 'Table of Contents' with links to various sections, including 'Abstract', 'Committee Membership', 'Foreword', 'Overview of Changes', and several chapters. The main content area has a blue header with the document title: 'CLSI M45-ED3:2016 Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria, 3rd Edition'. Below the header is a search box labeled 'Search within this Document' and a navigation bar with '< Previous | Next >' links. The main text describes the guideline's purpose and provides a list of authors.

CLSI M45 ED3:2016
Go To: [Top](#) | [Bottom](#)

Table of Contents -

- [Abstract](#)
- [Committee Membership](#)
- [Foreword](#)
- [Overview of Changes](#)
- [Chapter 1: Introduction](#)
- [Chapter 2: Indications for Performing Susceptibility Tests](#)
- [Chapter 3: Methods for Antimicrobial Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria](#)
- [Chapter 4: Quality System Essentials for Antimicrobial Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria](#)
- [Information and Interpretive Criteria for Susceptibility Testing](#)

Related Information +

CLSI M45-ED3:2016 Methods for Antimicrobial Dilution and Disk Susceptibility Testing of Infrequently Isolated or Fastidious Bacteria, 3rd Edition

Search within this Document

< Previous | Next >


This guideline informs clinical, public health, and research laboratories on susceptibility testing of infrequently isolated or fastidious bacteria that are not included in CLSI documents M02, M07, or M100. Antimicrobial agent selection, test interpretation, and quality control are addressed.

A guideline for global application developed through the Clinical and Laboratory Standards Institute consensus process.

M45, 3rd ed. August 2016 Replaces M45-A2

Janet A. Hindler, MCLS, MT(ASCP)	Romney M. Humphries, PhD, D(ABMM)
Sandra S. Richter, MD, D(ABMM)	James H. Jorgensen, PhD
Kathy Bernard, MSc, ARM(CCM)	Scott B. Killian, BA
Sonya Bodeis-Jones, BS	Peggy Kohner, BS, MT(ASCP)
Mariana Castanheira, PhD	Erika Matuschek, PhD
Diane M. Citron, BS, M(ASCP)	Patrick McDermott, PhD
Marc R. Couturier, PhD, D(ABMM)	Samir Patel, PhD, FCCM

IF FED UP WITH THE INTERNET



382 pages
± 30



Three General Comments



(NON-)FASTIDIOUS GROUPINGS

- Group A Primary test and report
- Group B Optional primary test, report selectively
- Group C Supplemental report selectively
- Group U Supplemental for urine only

TABLES 1

Table 1A
Suggested Nonfastidious Groupings
M02 and M07

Table 1A. Suggested Groupings of Antimicrobial Agents Approved by the US Food and Drug Administration for Clinical Use That Should Be Considered for Testing and Reporting on Nonfastidious Organisms by Microbiology Laboratories in the United States

	<i>Enterobacteriaceae</i>	<i>Pseudomonas aeruginosa</i>	<i>Staphylococcus</i> spp.	<i>Enterococcus</i> spp. ^m
--	---------------------------	-------------------------------	----------------------------	---------------------------------------

Table 1B
Suggested Fastidious Groupings
M02 and M07

Table 1B. Suggested Groupings of Antimicrobial Agents Approved by the US Food and Drug Administration for Clinical Use That Should Be Considered for Testing and Reporting on Fastidious Organisms by Microbiology Laboratories in the United States

ST	<i>Haemophilus influenzae</i> ^d and <i>Haemophilus parainfluenzae</i>	<i>Neisseria gonorrhoeae</i> ^l	<i>Streptococcus pneumoniae</i> ^l	<i>Streptococcus</i> spp. β-Hemolytic Group ^p	<i>Streptococcus</i> spp. Viridans Group ^p
----	---	---	--	---	--

Table 1C
Suggested Anaerobe Groupings
M11

Table 1C. Suggested Groupings of Antimicrobial Agents Approved by the US Food and Drug Administration for Clinical Use That Should Be Considered for Testing and Reporting on Anaerobic Organisms by Microbiology Laboratories in the United States

	Gram-Negative Anaerobes	Gram-Positive Anaerobes ^a
--	-------------------------	--------------------------------------

CRITERIA FOR INCLUSION

Agents of proven efficacy
Acceptable *in vitro* test performance

CRITERIA FOR ASSIGNMENT

Clinical efficacy
Prevalence of resistance
Minimizing emergence of resistance
FDA clinical indications for use
Current consensus recommendations for first-choice or alternative drugs
Co\$t


TABLE 1 GROUPINGS

- Tier 1 Antimicrobial agents that are appropriate for routine, primary testing and reporting
- Tier 2 Antimicrobial agents that are appropriate for routine, primary testing but may be reported following cascade reporting rules established at each institution



TABLE 1 GROUPINGS

- Tier 3 Antimicrobial agents that are appropriate for routine, primary testing in institutions that serve patients at high risk for MDROs but should only be reported following cascade reporting rules established at each institution



- Tier 4 Antimicrobial agents that may warrant testing and reporting by clinician request if antimicrobial agents in other tiers are not optimal because of various factors

REPORTING

- **Selective** Based on defined criteria unrelated to susceptibility testing data


Organism ID

Clinical setting

Site of infection

Patient demographics

- **Cascade** Based on overall antimicrobial susceptibility profile of isolate



CLSI M100-Ed33, 2023

THEY NOW PERFECTLY MATCH UP

Table 1B-2
Acinetobacter spp.
CLSI M02 and CLSI M07

Table 1B-2. *Acinetobacter* spp.

Tier 1: Antimicrobial agents that are appropriate for routine, primary testing and reporting	Tier 2: Antimicrobial agents that are appropriate for routine, primary testing but may be reported following cascade reporting rules established at each institution	Tier 3: Antimicrobial agents that are appropriate for routine, primary testing in institutions that serve patients at high risk for MDROs but should only be reported following cascade reporting rules established at each institution	Tier 4: Antimicrobial agents that may warrant testing and reporting by clinician request if antimicrobial agents in other tiers are not optimal because of various factors
--	--	---	--

Table 2B-2
Acinetobacter spp.
CLSI M02 and CLSI M07

Table 2B-2. Zone Diameter and MIC Breakpoints for *Acinetobacter* spp.

<p>Testing Conditions</p> <p>Medium: Disk diffusion: MHA Broth dilution: CAMHB; iron-depleted CAMHB for cefiderocol (see Appendix H)¹ Agar dilution: MHA</p> <p>Inoculum: Broth culture method or colony suspension, equivalent to a 0.5 McFarland standard; positive blood culture broth for select antimicrobial agents with disk diffusion (see general comment [3])</p> <p>Incubation: 35°C ± 2°C; ambient air; 20–24 hours, all methods</p>	<p>Routine QC Recommendations (see Tables 4A-1 and 5A-1 for acceptable QC ranges)</p> <p><i>Escherichia coli</i> ATCC^{®a} 25922 (for tetracyclines and trimethoprim-sulfamethoxazole)</p> <p><i>Pseudomonas aeruginosa</i> ATCC[®] 27853</p> <p>Refer to Tables 4A-2 and 5A-2 to select strains for routine QC of β-lactam combination agents.</p> <p>When a commercial test system is used for susceptibility testing, refer to the manufacturer's instructions for QC test recommendations and QC ranges.</p>
--	--

TABLE 2B-2

Table 2B-2. *Acinetobacter* spp. (Continued)

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, µg/mL			Comments
		S	I	R	S	I	R	
PENICILLINS								
Piperacillin*	100 µg	≥ 21	18–20	≤ 17	≤ 16	32–64	≥ 128	
β-LACTAM COMBINATION AGENTS								
(4) Organisms that test susceptible to the β-lactam agent alone are also considered susceptible to the β-lactam combination agent. However, organisms that test susceptible to the β-lactam combination agent cannot be assumed to be susceptible to the β-lactam agent alone. Similarly, organisms that test intermediate or resistant to the β-lactam agent alone may be susceptible to the β-lactam combination agent.								
Ampicillin-sulbactam	10/10 µg	≥ 15	12–14	≤ 11	≤ 8/4	16/8	≥ 32/16	
Piperacillin-tazobactam	100/10 µg	≥ 21	18–20	≤ 17	≤ 16/4	32/4–64/4	≥ 128/4	
Sulbactam-durlobactam	10/10 µg	≥ 17	14–16	≤ 13	≤ 4/4	8/4	≥ 16/4	
Ticarcillin-clavulanate*	75/10 µg	≥ 20	15–19	≤ 14	≤ 16/2	32/2–64/2	≥ 128/2	
CEPHEMS (PARENTERAL) (Including cephalosporins I, II, III, and IV. Please refer to Glossary I.)								
Ceftazidime	30 µg	≥ 18	15–17	≤ 14	≤ 8	16	≥ 32	
Cefepime	30 µg	≥ 18	15–17	≤ 14	≤ 8	16	≥ 32	
Cefotaxime	30 µg	≥ 23	15–22	≤ 14	≤ 8	16–32	≥ 64	
Ceftriaxone	30 µg	≥ 21	14–20	≤ 13	≤ 8	16–32	≥ 64	

- 1,2 A-1 Enterobacterales (excluding *Salmonella/Shigella*)
- 1,2 A-2 *Salmonella* and *Shigella* spp.
- 1,2 B-1 *Pseudomonas aeruginosa*
- 1,2 B-2 *Acinetobacter* spp.
- 1,2 B-3 *Burkholderia cepacia* complex
- 1,2 B-4 *Stenotrophomonas maltophilia*
- 1,2 B-5 Other Non-Enterobacterales
- 1,2 C *Staphylococcus* spp.
- 1,2 D *Enterococcus* spp.
- 1,2 E *Haemophilus influenzae* and *Haemophilus parainfluenzae*
- 1,2 F *Neisseria gonorrhoeae*
- 1,2 G *Streptococcus pneumoniae*
- 1,2 H-1 *Streptococcus* spp. β -Hemolytic Group
- 1,2 H-2 *Streptococcus* spp. Viridans Group
- 1,2 I *Neisseria meningitidis* (now has Table 1I)
- 1,2 J **Anaerobes** (combined Gram-positive and Gram-negative)

CLARIFICATIONS...or maybe not

- Daptomycin not routinely reported on organisms isolated from lower respiratory tract
- *Streptococcus agalactiae* intrapartum guidelines
- Susceptible isolates that may develop resistance after initiation of therapy

Formerly “within 3 to 4 days”

Now “within a few days”

BEFORE

Table 2D. *Enterococcus* spp. (Continued)

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, µg/mL				Comments
		S	I	R	S	SDD	I	R	
PENICILLINS									
Penicillin	10 units	≥15	-	≤14	≤8	-	-	≥16	<p>(7) The results of ampicillin susceptibility tests should be used to predict the activity of amoxicillin. Ampicillin results may be used to predict susceptibility to amoxicillin-clavulanate, ampicillin-sulbactam, and piperacillin-tazobactam among non-β-lactamase-producing enterococci. Ampicillin susceptibility can be used to predict imipenem susceptibility, providing the species is confirmed to be <i>E. faecalis</i>.</p> <p>(8) Enterococci susceptible to penicillin are predictably susceptible to ampicillin, amoxicillin, ampicillin-sulbactam, amoxicillin-clavulanate, and piperacillin-tazobactam for non-β-lactamase-producing enterococci. However, enterococci susceptible to ampicillin cannot be assumed to be susceptible to penicillin. If penicillin results are needed, testing of penicillin is required.</p> <p>(9) Rx: Combination therapy with high-dosage parenteral ampicillin, amoxicillin, penicillin, or vancomycin (for susceptible strains only), plus an aminoglycoside, is usually indicated for serious enterococcal infections, such as endocarditis, unless high-level resistance to both gentamicin and streptomycin is documented; such combinations are predicted to result in synergistic killing of enterococci.</p> <p>(10) Breakpoints are based on an ampicillin dosage regimen of 2 g parenterally administered every 4-6 h or an amoxicillin dosage regimen of 1-2 g parenterally administered every 6 h.</p> <p>(11) Breakpoints when oral ampicillin is used for therapy of uncomplicated UTIs only are based on an ampicillin dosage regimen of 500 mg orally administered every 6 h or amoxicillin dosage regimen of 250 mg orally administered every 8 h or 500 mg every 12 h.</p>
Ampicillin	10 µg	≥17	-	≤16	≤8	-	-	≥16	

Table 2D
Enterococcus spp.
M02 and M07

AFTER

Table 2D
Enterococcus spp.
 CLSI M02 and CLSI M07

Table 2D. *Enterococcus* spp. (Continued)

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, µg/mL				Comments
		S	I	R	S	SDD	I	R	
PENICILLINS (Continued)									
Penicillin	10 units	≥ 15	–	≤ 14	≤ 8	–	–	≥ 16	(10) Penicillin or ampicillin resistance among enterococci due to β-lactamase production has been reported very rarely. Penicillin or ampicillin resistance due to β-lactamase production is not reliably detected with routine disk or dilution methods but is detected using a direct, nitrocefin-based β-lactamase test. Because of the rarity of β-lactamase–positive enterococci, this test does not need to be performed routinely but can be used in selected cases. A positive β-lactamase test predicts resistance to penicillin as well as amino- and ureidopenicillins (see Glossary I).
Ampicillin	10 µg	≥ 17	–	≤ 16	≤ 8	–	–	≥ 16	

BEFORE

Table 3E-2. Enterobacterales (Continued)

Antimicrobial Agent	Disk Content	Read Times, hours	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm				Comments
			S	SDD	I	R	
PENICILLINS							
Ampicillin	10 µg	8-10	≥16	-	12-15	≤11	(4) Results of ampicillin testing can be used to predict results for amoxicillin. (5) Breakpoints are based on an ampicillin dosage regimen of 2 g parenterally administered every 4-6 h or an amoxicillin dosage regimen of 1-2 g parenterally administered every 6 h.
		16-18	≥17	-	14-16	≤13	
CEPHEMS (PARENTERAL) (Including cephalosporins I, II, III, and IV. Please refer to Glossary I.)							
Ceftriaxone	30 µg	8-10	≥23	-	20-22	≤19	(6) Breakpoints are based on a dosage regimen of 1 g administered every 24 h.
		16-18	≥23	-	20-22	≤19	
Ceftazidime	30 µg	8-10	≥21	-	18-20	≤17	(7) Breakpoints are based on a dosage regimen of 1 g administered every 8 h.
		16-18	≥21	-	18-20	≤17	
MONOBACTAMS							
Aztreonam	30 µg	8-10	≥21	-	18-20	≤17	(8) Breakpoints are based on a dosage regimen of 1 g administered every 8 h.
		16-18	≥21	-	18-20	≤17	

Table 3E-2
Zone Diameter Disk Diffusion Breakpoints for
Enterobacterales Direct From Blood Culture

AFTER

Table 3F-2
Zone Diameter Disk Diffusion Breakpoints for Enterobacterales Direct From Blood Culture

Table 3F-2. Zone Diameter Disk Diffusion Breakpoints for Enterobacterales Direct From Blood Culture

General Comments

- (1) Organism identification must be known before interpreting and reporting results. Fluoroquinolone breakpoints do not apply to *Salmonella* spp. Aztreonam, ceftazidime, and tobramycin breakpoints do not apply to *Salmonella* or *Shigella* spp.
- (2) For additional testing and reporting recommendations, refer to Tables 2A-1 and 2A-2.

NOTE: Information in boldface type is new or modified since the previous edition.

Antimicrobial Agent	Disk Content	Read Times, hours	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Comments
			S	I	R	
PENICILLINS						
Ampicillin	10 µg	8–10	≥ 16	12–15	≤ 11	(3) Results of ampicillin testing can be used to predict results for amoxicillin.
		16–18	≥ 17	14–16	≤ 13	
CEPHEMS (PARENTERAL) (Including cephalosporins I, II, III, and IV. Please refer to Glossary I.)						
Ceftriaxone	30 µg	8–10	≥ 23	20–22	≤ 19	
		16–18	≥ 23	20–22	≤ 19	
Ceftazidime	30 µg	8–10	≥ 21	18–20	≤ 17	
		16–18	≥ 21	18–20	≤ 17	
MONOBACTAMS						
Aztreonam	30 µg	8–10	≥ 21	18–20	≤ 17	
		16–18	≥ 21	18–20	≤ 17	

WHERE DID THEY GO?

Introduction to Table 2 Dosages.
Dosage Regimens Used to Establish Susceptible or Susceptible-Dose
Dependent Breakpoints

Introduction to Table 2 Dosages. Antimicrobial Agent Dosage Regimens Used to Establish Susceptible or Susceptible-Dose Dependent Breakpoints

The evolving science of pharmacokinetics/pharmacodynamics has become increasingly important in recent years in determining MIC breakpoints. CLSI susceptible or susceptible-dose dependent breakpoints **added or revised since 2010** have been based on a specific dosage regimen(s); these dosage regimens are listed in the table below. Proper application of the breakpoints necessitates drug exposure at the site of infection that corresponds to or exceeds the expected systemic drug exposure at the dose listed in adult patients with normal renal function. This information should be shared with pharmacists, infectious diseases staff, and others making dosing recommendations for the institution.

CLSI guidance for establishing or revising breakpoints is available in CLSI M23.¹ Rationale documents that provide the scientific reasoning behind the subcommittee's decisions for some breakpoints, along with documentation of the standardized data and methods used to determine breakpoints, can be found on the CLSI website.²

NOTE 1: If both a susceptible and a susceptible-dose dependent dosage regimen were used, they are designated by "S" or "SDD" preceding the dosage regimen. Otherwise, it should be assumed that the dosage regimen applies to the susceptible breakpoint.

NOTE 2: Unless otherwise noted, refer to the approved prescribing information for the infusion duration used to set breakpoints for IV antibiotics (eg, 0.5 hours for most β -lactams, 1–1.5 hours for fluoroquinolones).

NOTE 3: Dosage regimens also include the frequency of administration designated by the abbreviation "q." For example, the amikacin susceptible breakpoint for Enterobacterales was based on a dosage regimen of 15 mg/kg IV q 24 h, which corresponds to 15 mg/kg IV administered every 24 hours.

SUSCEPTIBLE AND -DOSE DEPENDENT

Table 2A-1
Enterobacterales (excluding *Salmonella/Shigella*)
CLSI M02 and CLSI M07

Table 2A-1. Enterobacterales (excluding *Salmonella/Shigella*) (Continued)

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm				Interpretive Categories and MIC Breakpoints, µg/mL				Comments
		S	SDD	I	R	S	SDD	I	R	
CEPHEMS (PARENTERAL) (Including cephalosporins I, II, III, and IV. Please refer to Glossary I.) (Continued)										
Cefepime	30 µg	≥ 25	19–24	–	≤ 18	≤ 2	4–8	–	≥ 16	(18) Cefepime S/SDD results should be suppressed or edited and reported as resistant for isolates that demonstrate carbapenemase production (see Appendix G, Table G3).

Table 2 Dosages. Dosage Regimens Used to Establish Susceptible or Susceptible-Dose Dependent Breakpoints
CLSI M02 and CLSI M07

Table 2 Dosages. Antimicrobial Agent Dosage Regimens Used to Establish Susceptible or Susceptible-Dose Dependent Breakpoints

Antimicrobial Agent	Dosage Regimen Used to Establish S or SDD Breakpoint
Table 2A-1. Enterobacterales (excluding <i>Salmonella/Shigella</i>)	
Amikacin	15 mg/kg IV q 24 h
Ampicillin (ampicillin test results predict results for amoxicillin)	Ampicillin: 2 g IV q 4–6 h or Amoxicillin: 1–2 g IV q 6 h
Ampicillin (ampicillin test results predict results for amoxicillin; <i>Escherichia coli</i> and <i>Proteus mirabilis</i> for uncomplicated UTIs only)	Ampicillin: 500 mg PO q 6 h or Amoxicillin: 250 mg PO q 8 h or 500 mg PO q 12 h
Amoxicillin-clavulanate (oral amoxicillin-clavulanate for uncomplicated UTIs or when completing therapy for systemic infection only)	1.2 g (1 g amoxicillin + 0.2 g clavulanate) IV q 6 h 500/125 mg PO q 8 h or 875/125 mg PO q 12 h
Ampicillin-sulbactam	3 g IV (2 g ampicillin + 1 g sulbactam) q 6 h
Aztreonam	1 g IV q 8 h
Cefazolin (<i>E. coli</i> , <i>Klebsiella pneumoniae</i> , and <i>P. mirabilis</i> for infections other than uncomplicated UTIs only)	2 g IV q 8 h
Cefazolin (<i>E. coli</i> , <i>K. pneumoniae</i> , and <i>P. mirabilis</i> for uncomplicated UTIs only)	1 g IV q 12 h
Ceftaroline	600 mg IV q 12 h
Cefepime	S: 1 g IV q 8 h or 2 g IV q 12 h SDD: 2 g IV q 8 h over 3 h

HERE'S OUR *Enterococcus*



Table 2D. *Enterococcus* spp.

Ampicillin (ampicillin test results predict results for amoxicillin; oral ampicillin or amoxicillin used for uncomplicated UTIs only)	Ampicillin: 2 g IV q 4–6 h or 500 mg PO q 6 h Amoxicillin: 1–2 g IV q 6 h or 250 mg PO q 8 h or 500 mg PO q 12 h
Dalbavancin (vancomycin-susceptible <i>E. faecalis</i> only)	1500 mg IV once or 1000 mg IV once followed one week later by 500 mg IV once
Daptomycin (<i>E. faecium</i> only)	SDD: 8–12 mg/kg IV q 24 h
Daptomycin (<i>Enterococcus</i> spp. other than <i>E. faecium</i>)	6 mg/kg IV q 24 h
Oritavancin (vancomycin-susceptible <i>E. faecalis</i> only)	1200 mg IV once
Tedizolid (<i>E. faecalis</i> only)	200 mg IV/PO q 24 h
Telavancin (vancomycin-susceptible <i>E. faecalis</i> only)	10 mg/kg IV q 24 h



Table 2E. *Haemophilus influenzae* and *Haemophilus parainfluenzae*

Amoxicillin-clavulanate	500/125 mg PO q 8 h or 875/125 mg PO q 12 h
Ampicillin (meningitis)	2 g IV q 4 h
Ampicillin-sulbactam	3 g (2 g ampicillin + 1 g sulbactam) IV q 6 h
Ceftaroline (<i>H. influenzae</i> only)	600 mg IV q 12 h
Ceftolozane-tazobactam (<i>H. influenzae</i> only)	3 g (2 g ceftolozane + 1 g tazobactam) IV q 8 h

Table 2 Dosages. Dosage Regimens Used to Establish Susceptible or Susceptible-Dose Dependent Breakpoints
CLSI M02 and CLSI M07



Old Business



Table 3G-1
Test for Detecting Methicillin (Oxacillin) Resistance in *Staphylococcus aureus* and *Staphylococcus lugdunensis*

Table 3G-1. Test for Detecting Methicillin (Oxacillin) Resistance in *Staphylococcus aureus*^a and *Staphylococcus lugdunensis*

Test	Detecting <i>mecA</i> -Mediated Resistance Using Cefoxitin ^b		Detecting <i>mecA</i> -Mediated Resistance Using Oxacillin	Detecting <i>mecA</i> -mediated Resistance Using Oxacillin Salt Agar for <i>S. aureus</i> Only
Test method	Disk diffusion	Broth microdilution	Broth microdilution and agar dilution	Agar dilution for <i>S. aureus</i>
Medium	MHA	CAMHB	CAMHB with 2% NaCl (broth microdilution) MHA with 2% NaCl (agar dilution)	MHA with 4% NaCl
Antimicrobial concentration	30-µg cefoxitin disk	4 µg/mL cefoxitin	2 µg/mL oxacillin	6 µg/mL oxacillin
Inoculum	Standard disk diffusion procedure	Standard broth microdilution procedure	Standard broth microdilution procedure or standard agar dilution procedure	Colony suspension to obtain 0.5 McFarland turbidity Using a 1-µL loop that was dipped in the suspension, spot an area 10-15 mm in diameter. Alternatively, using a swab dipped

Table 3H
Oxacillin Salt Agar Test for Methicillin (Oxacillin) Resistance in *Staphylococcus aureus*

Table 3H. Oxacillin Salt Agar Test for Detecting Methicillin (Oxacillin) Resistance in *Staphylococcus aureus*^a

Test	Oxacillin Salt Agar
Test method	Agar dilution
Medium	MHA with 4% NaCl
Antimicrobial concentration	6 µg/mL oxacillin
Inoculum	Colony suspension to obtain 0.5 McFarland turbidity Using a 1-µL loop that was dipped in the suspension, spot an area 10–15 mm in diameter. Alternatively, using a swab dipped in the suspension and the excess liquid expressed, spot a similar area or streak an entire quadrant.

THE STREAK COMES TO AN END

Table 2C
Staphylococcus spp.
CLSI M02 and CLSI M07


Table 2C. *Staphylococcus* spp. (Continued)

Methods or Targets for Detection of Methicillin (Oxacillin)-Resistant <i>Staphylococcus</i> spp.							
Organism	Disk Diffusion		MIC		<i>mecA</i>	PBP2a	Oxacillin Salt Agar
	Cefoxitin	Oxacillin	Cefoxitin	Oxacillin			
<i>S. aureus</i>	Yes (16–18 h)	No	Yes (16–20 h)	Yes (24 h)	Yes	Yes	Yes (24 h)
<i>S. lugdunensis</i>	Yes (16–18 h)	No	Yes (16–20 h)	Yes (24 h)	Yes	Yes	No
<i>S. epidermidis</i>	Yes (24 h)	Yes (16–18 h)	No	Yes (24 h)	Yes	Yes	No
<i>S. pseudintermedius</i>	No	Yes (16–18 h)	No	Yes (24 h)	Yes	Yes	No
<i>S. schleiferi</i>	No	Yes (16–18 h)	No	Yes (24 h)	Yes	Yes	No
<i>Staphylococcus</i> spp. (not listed above or not identified to the species level)	Yes, with exceptions ^a (24 h)	No	No	Yes (24 h)	Yes	Yes	No

Abbreviations: h, hour(s); MIC, minimal inhibitory concentration; PBP2a, penicillin-binding protein 2a.

^a The cefoxitin disk diffusion test may not perform reliably for all species (eg, *S. haemolyticus*) that fall into the category of “*Staphylococcus* spp. (not listed above or not identified to the species level).”⁶

- Grouped by method, rather than test (data are the same)
- Added *mecA* and PBP2a determinations



AMERICAN
SOCIETY FOR
MICROBIOLOGY

Journal of
Clinical Microbiology®

COMMENTARY



Use of Rapid Diagnostics To Manage Pediatric Bloodstream Infections? You Bet Your ASP!

Mark D. Gonzalez,^a Melanie L. Yarbrough^b

^aDepartment of Pathology and Laboratory Services, Children's Healthcare of Atlanta, Atlanta, Georgia, USA

^bDepartment of Pathology & Immunology, Washington University School of Medicine, St. Louis, Missouri, USA

POSITIVE BLOOD CULTURE BOTTLE

JOURNAL OF CLINICAL MICROBIOLOGY, Mar. 1979, p. 347-350
0095-1137/79/03-0347/04\$02.00/0

Vol. 9, No. 3


Standardization of Direct Susceptibility Test for Blood Cultures

DALE FAY†* AND JEAN E. OLDFATHER

Riverside Methodist Hospital, Columbus, Ohio 43214

Received for publication 17 December 1978

Insufficient data are available to establish the reliability of direct disk diffusion susceptibility tests performed utilizing positive blood culture broth as inoculum. When *Staphylococcus aureus* ATCC 25923, *Escherichia coli* ATCC 25922, and *Pseudomonas aeruginosa* ATCC 27853 were used, 0.03 ml of turbid overnight blood culture broth was found to produce zone diameters closely approximating the size of diameters obtained by a standardized method. Results of direct (0.03 ml of inoculum) and standardized susceptibility tests were then compared for 116 positive blood cultures (1,069 individual disk comparisons). There were 1,011 test agreements (94.6%). There were also 48 (4.5%) minor discrepancies (change between sensitive and intermediate or between intermediate and resistant) and 10 (0.9%) major discrepancies (change between sensitive and resistant). The major discrepancies were randomly distributed among several organisms and antibiotics. Discrepancies occurred most frequently in the more clinically acceptable direction; i.e., in 79.3% the direct test indicted greater resistance than the standardized test. These data establish that 0.03 ml of turbid overnight blood culture broth produces results which compare closely to those obtained with standard methods, and in practice yield direct susceptibility results with a clinically acceptable level of reliability.



TWO DROPS

TABLE 2. Organisms included in the clinical comparison of the direct and standardized susceptibility tests

Organism	No. of strains tested	Discrepancies		Agreements
		Major	Minor	
<i>E. coli</i>	46	2	22	390
<i>Klebsiella</i>	16	2	12	130
<i>Proteus mirabilis</i>	8	2	3	67
<i>Providencia stuartii</i>	1	0	1	8
<i>Citrobacter diversus</i>	1	0	0	9
<i>Citrobacter freundii</i>	1	0	1	8
<i>Enterobacter aerogenes</i>	3	0	1	26
<i>Enterobacter cloacae</i>	3	2	1	24
<i>Enterobacter agglomerans</i>	1	0	0	9
<i>Serratia marcescens</i>	3	0	0	27
<i>P. aeruginosa</i>	4	0	0	36
<i>Pseudomonas</i> species	2	0	2	16
<i>Bordetella pertussis</i>	1	0	0	9
<i>Acinetobacter calcoaceticus</i>	1	0	0	9
<i>S. aureus</i>	12	0	0	120
<i>Staphylococcus epidermidis</i>	8	2	4	74
Enterococcus	3	0	1	29
Group D <i>Streptococcus</i> (not Enterococcus)	1	0	0	10
Viridans <i>Streptococcus</i>	1	0	0	10

TABLE 3. Distribution of discrepancies between direct and standardized susceptibility tests by antibiotic

Antibiotic	No. of comparisons	Discrepancies		
		Total	Major	Minor
Ampicillin	116	4 (3.8) ^a	1	3
Carbenicillin	91	4 (4.3)	0	4
Cephalothin	116	16 (13.8)	2	14
Chloramphenicol	116	6 (5.2)	3	3
Clindamycin	25	0	0	0
Colistin	91	6 (6.6)	2	4
Erythromycin	25	0	0	0
Gentamicin	116	0	0	0
Kanamycin	116	1 (0.8)	0	1
Methicillin	25	1 (4.0)	0	1
Penicillin	25	3 (12.0)	1	2
Streptomycin	91	9 (9.9)	0	9
Tetracycline	116	8 (6.9)	1	7

Major (0.9%): shift between sensitive and resistant

Minor (4.5%): shift between sensitive and intermediate
shift between intermediate and resistant

FUTHER ENHANCEMENT

JOURNAL OF CLINICAL MICROBIOLOGY, Sept. 1984, p. 473-477
0095-1137/84/090473-05\$02.00/0
Copyright © 1984, American Society for Microbiology

Vol. 20, No. 3

Rapid Antimicrobial Susceptibility Testing of Isolates from Blood Cultures by Direct Inoculation and Early Reading of Disk Diffusion Tests

MARIE B. COYLE,^{1,2*} LEE ANNE MCGONAGLE,³ JAMES J. PLORDE,⁴ CARLA R. CLAUSEN,⁵ AND FRITZ D. SCHOENKNECHT³

Clinical Microbiology Division, University of Washington,¹ and University Hospital,³ Seattle, Washington 98195; Harborview Medical Center, Seattle, Washington 98104;² Seattle Veterans Administration Medical Center, Seattle, Washington 98108;⁴ and Childrens Orthopedic Hospital and Medical Center, Seattle, Washington 98105⁵

Received 13 January 1984/Accepted 24 May 1984

Disk diffusion tests, inoculated directly from positive blood cultures, were evaluated for accuracy of reading zone diameters after 4- and 6-h and overnight incubation. In comparisons with results from standard disk diffusion tests, the 4-h results were in agreement for 83% of tests with gram-positive organisms and 64% of tests with gram-negative organisms. When minor discrepancies were ignored, the 4-h readings were in agreement for 98% of the tests with gram-positive organisms and 95% of the tests with gram-negative organisms. After 6 h of incubation, 91% of the tests with gram-positive organisms and 86% of the tests with gram-negative organisms agreed with standard results. The agreement was 99% for tests with both gram-positive and gram-negative organisms when minor discrepancies were excluded. Very major discrepancies occurred in two tests (0.1%) with gram-positive organisms and were not observed in tests with gram-negative organisms. The frequencies of major discrepancies were 3.5% after 4 h, 0.6% after 6 h, and 0.7% after overnight incubation. Ampicillin and cephalothin tests with *Escherichia coli* and *Klebsiella* spp. accounted for 81% of the major discrepancies in tests with gram-negative organisms. Oxacillin tests accounted for more than half of the major discrepancies in tests with staphylococci. The results of this study, which did not include the newer antibiotics, indicate that direct susceptibility tests from blood cultures read after 6 h of incubation are more reliable than 4-h results and produce less than 1% major errors in comparisons with standard susceptibility tests.



READING 'EM EARLY

TABLE 1. Percentage of isolates with direct tests read after 4 or 6 h

Blood culture isolate	No. of isolates	% Read after:	
		4 h	6 h ^a
Gram positive			
<i>S. aureus</i>	60	22	63
Coagulase-negative staphylococci	87	3	21
Beta-hemolytic streptococci	30	37	87
Enterococci	21	19	52
Pneumococci	21	10	38
Viridans streptococci	14	0	0
Total for gram positive	233	14	44
Gram negative			
<i>E. coli</i>	84	52	85
<i>Klebsiella spp.</i>	38	40	76
<i>Enterobacter spp.</i>	12	42	92
<i>P. aeruginosa</i>	11	0	64
Others ^b	25	36	60
Total for gram negative	170	43	78

TABLE 2. Discrepancies from direct tests compared with standardized tests

Isolate type, time incubated	No. of tests	No. of discrepancies			Overall agreement (%)
		Very major (%)	Major (%)	Minor (%)	
Gram positive					
4 h	216	1 (0.5)	3 (1.4)	32 (14.8)	83.3
6 h	494	0	3 (0.6)	39 (7.9)	91.4
Overnight	1,307	1 (0.07)	8 (0.6)	65 (5.0)	94.3
Gram negative					
4 h	361	0	17 (4.7)	114 (31.6)	63.7
6 h	438	0	3 (0.7)	59 (13.5)	85.8
Overnight	762	0	6 (0.8)	73 (9.6)	89.6

four



AMERICAN
SOCIETY FOR
MICROBIOLOGY

Journal of
Clinical Microbiology®

BACTERIOLOGY



Direct-from-Blood-Culture Disk Diffusion To Determine Antimicrobial Susceptibility of Gram-Negative Bacteria: Preliminary Report from the Clinical and Laboratory Standards Institute Methods Development and Standardization Working Group

- Resistance in GNR can be multi-factorial; full phenotypic approach may be desirable
- Little standardization; very few laboratories report
- 1 carbapenem-resistant *Acinetobacter baumannii*

J Clin Microbiol. 56:e01678-17; 2018

RESULTS

TABLE 3 Resolved performance of direct-from-blood-culture disk diffusion method at 18 h, by antibiotic

Drug	No. of isolates		% CA	No. (%) of:		
	S	R		VME	ME	mE
Amikacin	45	13	96.7	0 (0)	0 (0)	2 (3.3)
Amoxicillin-clavulanate	9	17	88.9	0 (0)	1 (11.1)	2 (7.4)
Ampicillin ←	6	9	93.3	0 (0)	0 (0)	1 (6.7)
Aztreonam ←	21	28	94.3	0 (0)	0 (0)	3 (5.7)
Cefazolin	5	18	73.1	0 (0)	2 (40.0)	5 (19.2)
Cefepime ←	41	17	91.7	0 (0)	0 (0)	5 (8.3)
Cefoxitin	10	15	85.2	0 (0)	1 (10.0)	3 (11.1)
Ceftazidime ←	25	31	89.8	0 (0)	0 (0)	6 (10.2)
Ceftriaxone ←	16	29	87.5	0 (0)	2 (12.5)	4 (8.3)
Ciprofloxacin ←	26	27	96.6	0 (0)	0 (0)	1 (1.7)
Ertapenem	22	12	83.3	0 (0)	0 (0)	7 (16.7)
Gentamicin	39	18	95.0	0 (0)	1 (2.6)	2 (3.3)
Imipenem	34	21	68.3	0 (0)	3 (8.8)	15 (25.0)
Levofloxacin	33	25	91.7	0 (0)	1 (3.0)	3 (5.0)
Meropenem ←	37	19	84.7	0 (0)	1 (2.7)	8 (13.6)
Minocycline	29	11	80.0	0 (0)	0 (0)	9 (20.0)
Piperacillin-tazobactam	23	30	83.3	0 (0)	0 (0)	10 (16.7)
Tigecycline	35	3	87.2	0 (0)	0 (0)	5 (12.8)
Tobramycin ←	39	17	93.2	0 (0)	0 (0)	4 (6.8)
Trimethoprim-sulfamethoxazole ←	17	30	95.8	0 (0)	0 (0)	2 (4.2)



RESULTS

TABLE 3 Resolved performance of cod-culture disk diffusion method at 6 h, 18 h, by antibiotic

Drug	% CA	No. (%) of:	
		VME	ME
Amikacin	96.7	0 (0)	0 (0)
Amoxicillin-clavulanate	88.9	0 (0)	1 (11.1)
Ampicillin	93.3	0 (0)	0 (0)
Aztreonam	94.3	0 (0)	0 (0)
Cefazolin	73.1	0 (0)	2 (40.0)
Cefepime	91.7	0 (0)	0 (0)
Cefoxitin	85.2	0 (0)	1 (10.0)
Ceftazidime	89.8	0 (0)	0 (0)
Ceftriaxone	87.5	0 (0)	2 (12.5)
Ciprofloxacin	96.6	0 (0)	0 (0)
Ertapenem	83.3	0 (0)	0 (0)
Gentamicin	95.0	0 (0)	1 (2.6)
Imipenem	68.3	0 (0)	3 (8.8)
Levofloxacin	91.7	0 (0)	1 (3.0)
Meropenem	84.7	0 (0)	1 (2.7)
Minocycline	80.0	0 (0)	0 (0)
Piperacillin-tazobactam	83.3	0 (0)	0 (0)
Tigecycline	87.2	0 (0)	0 (0)
Tobramycin	93.2	0 (0)	0 (0)
Trimethoprim-sulfamethoxazole	95.8	0 (0)	0 (0)

TABLE 5 Resolved performance of direct-from-blood-culture disk diffusion method at 6 h, by antibiotic

Drug	No. of isolates		% CA	No. (%) of:		
	S	R		VME	ME	mE
Amikacin	45	13	62.2	3 (23.1)	2 (4.4)	12 (26.7)
Amoxicillin-clavulanate	9	17	60.0	0 (0)	1 (11.1)	9 (36.0)
Ampicillin ←	6	9	69.2	0 (0)	1 (16.7)	3 (23.1)
Aztreonam ←	21	28	84.2	0 (0)	1 (4.8)	5 (13.2)
Cefazolin	5	18	66.7	1 (5.6)	2 (40.0)	6 (25.0)
Cefepime ←	41	17	75.6	0 (0)	4 (9.8)	6 (13.3)
Cefoxitin	10	15	68.0	0 (0)	1 (10.0)	7 (28.0)
Ceftazidime ←	25	31	65.9	0 (0)	4 (16.0)	11 (25.0)
Ceftriaxone ←	16	29	77.3	0 (0)	3 (18.8)	7 (15.9)
Ciprofloxacin ←	24	27	57.1	0 (0)	1 (4.2)	16 (39.0)
Ertapenem	22	12	73.7	0 (0)	2 (9.1)	8 (21.1)
Gentamicin	39	18	95.6	0 (0)	0	2 (4.4)
Imipenem	34	21	46.7	0 (0)	6 (17.6)	18 (40.0)
Levofloxacin	33	25	75.6	0 (0)	1 (3.0)	10 (22.2)
Meropenem ←	36	19	52.3	0 (0)	9 (25.0)	11 (25.6)
Minocycline	29	11	65.9	0 (0)	0	12 (29.3)
Piperacillin-tazobactam	22	30	64.4	2 (6.7)	4 (18.2)	11 (25.0)
Tigecycline	35	3	45.7	0 (0)	3 (8.6)	16 (45.7)
Tobramycin ←	39	17	95.6	0 (0)	0	2 (4.4)
Trimethoprim-sulfamethoxazole ←	17	30	86.4	1 (3.3)	2 (11.8)	3 (6.8)



Evaluation of the Performance of Direct Susceptibility Test by VITEK-2 from Positively Flagged Blood Culture Broth for Gram-Negative Bacilli

Kavipriya D¹ Suman Susan Prakash¹ Sarumathi Dhandapani¹ Deepashree Rajshekar¹
Apurba Sankar Sastry¹




Table 3 Performance of direct test compared with reference (colony) test for nonfermenters by VITEK-2 system

Nonfermenter (60)	Categorical agreement (%)	Categorical disagreement (%)				Essential agreement	
		Minor	Major	Very major	Total	Agreed	Disagreed
Ticarillin/clavulanic acid	57 (95.0%)	3 (5.0%)	0 (0.0)	0 (0.0)	3 (5.0%)	57 (95.0%)	3 (5.0%)
Piperacillin/tazobactam	59 (98.3%)	0 (0.0)	1 (1.7%)	0 (0.0)	1 (1.7%)	59 (98.3%)	1 (1.7%)
Ceftazidime	60 (100.0%)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	60 (100.0%)	0 (0.0)
Cefoperazone/sulbactam	59 (98.3%)	1 (1.7%)	0 (0.0)	0 (0.0)	1 (1.7%)	59 (98.3%)	1 (1.7%)
Cefepime	60 (100.0%)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	60 (100.0%)	0 (0.0)
Doripenem	60 (100.0%)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	60 (100.0%)	0 (0.0)
Imipenem	60 (100.0%)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	60 (100.0%)	0 (0.0)
Meropenem	60 (100.0%)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	60 (100.0%)	0 (0.0)
Amikacin	59 (98.3%)	0 (0.0)	0 (0.0)	1 (1.7%)	1 (1.7%)	59 (98.3%)	1 (1.7%)
Gentamicin	59 (98.3%)	1 (1.7%)	0 (0.0)	0 (0.0)	1 (1.7%)	59 (98.3%)	1 (1.7%)
Ciprofloxacin	59 (98.3%)	1 (1.7%)	0 (0.0)	0 (0.0)	1 (1.7%)	59 (98.3%)	1 (1.7%)
Levofloxacin	55 (91.7%)	5 (8.3%)	0 (0.0)	0 (0.0)	5 (8.3%)	55 (91.7%)	5 (8.3%)
Minocycline	54 (90.0%)	6 (10.0%)	0 (0.0)	0 (0.0)	6 (10.0%)	54 (90.0%)	6 (10.0%)
Tigecycline	57 (95.0%)	2 (3.3%)	1 (1.6%)	0 (0.0)	3 (5.0%)	57 (95.0%)	3 (5.0%)
Colistin	60 (100.0%)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	60 (100.0%)	0 (0.0)

[up to 15%
categorical
disagreement
(mode 1.7%) for
Enterobacterales]

Table 3F-1. Test for Performing Disk Diffusion Directly From Positive Blood Culture Broth


Test	Direct Disk Diffusion
Test method	Disk diffusion using positive blood culture broth
Organism group	Enterobacterales, <i>Pseudomonas aeruginosa</i> , and <i>Acinetobacter</i> spp.
Medium	MHA
Antimicrobial concentration	Standard disk contents for the antimicrobial agents are detailed in Table 3F-2 (Enterobacterales), Table 3F-3 (<i>P. aeruginosa</i>), and Table 3F-4 (<i>Acinetobacter</i> spp.).
Inoculum	Positive blood culture broth with gram-negative bacilli, used within 8 hours of flagging positive by the blood culture system
Test procedure	<ol style="list-style-type: none"> 1. Invert blood culture bottle 5–10 times to thoroughly mix. 2. Sterilize the top of the bottle with an alcohol wipe (allow to dry) and insert 20-gauge venting needle into the blood culture bottle. 3. Dispense 4 drops of blood culture broth onto an MHA plate. As a purity check, use an inoculated blood agar plate streaked for isolation. 4. Spread blood culture broth across the entire surface of the MHA plate using a sterile cotton swab. 5. Repeat this procedure by streaking twice more, rotating the plate approximately 60 degrees each time to ensure an even distribution of inoculum. 6. Leave the lid ajar for 3–5 minutes (ideally) but no more than 15 minutes. 7. Dispense antimicrobial disks onto the surface of the inoculated MHA plate. 8. Press each disk down to ensure complete contact with the agar surface. 9. Invert the plate and place in the incubator within 15 minutes of disks being applied.
Incubation conditions	35°C ± 2°C; ambient air
Incubation length	8–10 hours or 16–18 hours (refer to Tables 3F-2, 3F-3, and 3F-4 for antimicrobial agent-specific incubation lengths)
Results	<ol style="list-style-type: none"> 1. Examine the blood agar purity plate to ensure pure growth. 2. Examine the test plate to ensure confluent lawn of growth appropriate to read disk zone tests per CLSI M02.¹ 3. Measure the zone diameters according to routine disk diffusion recommendations in CLSI M02.¹ 4. Interpret results using the zone diameter breakpoints in Tables 3F-2, 3F-3, and 3F-4 if the gram-negative bacillus tested is confirmed to be an Enterobacterales, <i>P. aeruginosa</i>, or <i>Acinetobacter</i> spp., respectively. If species is identified as another organism, do not interpret or report results. 5. Report only the interpretive category and not the measured zone size.

4

35 ± 2

8-10 or

16-18



Daily or weekly QC; *E. coli* ATCC 25922, *P. aeruginosa* ATCC 27853, *E. coli* ATCC 35218 (NEW), others if necessary (NEW)

Acinetobacter spp. (TABLE 3F-4)

Agent	Concentration	Incubation time (hours)	Zone Diameters (mm)		
			S	I	R
ampicillin-sulbactam	10/10 µg	16-18	≥ 15	12-14	≤ 11
ceftazidime	30 µg	16-18	≥ 17	15-16	≤ 14
cefepime	30 µg	8-10	≥ 18	15-17	≤ 14
		16-18	≥ 18	15-17	≤ 14
ceftriaxone	30 µg	8-10	≥ 21	14-20	≤ 13
		16-18	≥ 20	13-19	≤ 12
meropenem	10 µg	8-10	≥ 18	15-17	≤ 14
		16-18	≥ 18	15-17	≤ 14
tobramycin	10 µg	8-10	≥ 15	13-14	≤ 12
		16-18	≥ 15	13-14	≤ 12
ciprofloxacin	5 µg	8-10	≥ 21	16-20	≤ 15
		16-18	≥ 21	16-20	≤ 15
trimethoprim-sulfamethoxazole	1.25/23.75 µg	8-10	≥ 16	11-15	≤ 10
		16-18	≥ 16	11-15	≤ 10

TABLE 3F-2 REVISIONS

Agent	Concentration	Incubation time (hours)	Zone Diameters (mm)		
			S	I	R
tobramycin	10 µg	8-10	≥ 17	13-16	≤ 12
		16-18	≥ 17	13-16	≤ 12

Intermediate increased


- * direct aztreonam, ceftazidime, tobramycin not for *Salmonella* or *Shigella* spp.
- * direct ciprofloxacin not for *Salmonella* spp.

Glucose fermenters
 Reduce nitrates to nitrites
 Non-spore-forming GNR
 Grows on routine media
 Facultative
 Oxidase-negative (except *Plesiomonas*)

TABLE 3F-3 NEWBIES

Agent	Concentration	Incubation time (hours)	Zone Diameters (mm)		
			S	I	R
cefepime	30 μ g	16-18	≥ 18	15-17	≤ 14
tobramycin	10 μ g	8-10	≥ 19	13-18	≤ 12
		16-18	≥ 19	13-18	≤ 12

* confirmatory cefepime MIC testing for zone diameters 15-17 mm



GETTING BUSIER EVERY YEAR

Table 3F-1
Test for Performing Disk Diffusion Directly From Positive Blood Culture Broth

Table 3F-1. (Continued)

Breakpoint Additions Since 2021 (Continued)	Antimicrobial Agent	Date of Addition (M100 Edition)	8–10 h	16–18 h
	<i>Pseudomonas aeruginosa</i>			
	Cefepime	February 2024 (M100-Ed34)		X
	Ceftazidime	February 2022 (M100-Ed32)		X
	Ciprofloxacin	February 2022 (M100-Ed32)	X	X
	Meropenem	February 2022 (M100-Ed32)		X
		March 2023 (M100-Ed33)	X	
	Tobramycin	February 2022 (M100-Ed32)	X	X
	<i>Acinetobacter</i> spp.			
	Ampicillin-sulbactam	February 2024 (M100-Ed34)		X
	Cefepime	February 2024 (M100-Ed34)	X	X
	Ceftazidime	February 2024 (M100-Ed34)		X
	Ceftriaxone	February 2024 (M100-Ed34)	X	X
	Ciprofloxacin	February 2024 (M100-Ed34)	X	X
	Meropenem	February 2024 (M100-Ed34)	X	X
	Tobramycin	February 2024 (M100-Ed34)	X	X
	Trimethoprim-sulfamethoxazole	February 2024 (M100-Ed34)	X	X
Breakpoint Revisions Since 2021	Enterobacterales			
	Tobramycin	February 2024 (M100-Ed34)	X	X
	<i>Pseudomonas aeruginosa</i>			
	Tobramycin	February 2024 (M100-Ed34)	X	X

Abbreviations: ATCC®, American Type Culture Collection; MHA, Mueller-Hinton agar; QC, quality control.

CLSI M100-Ed34, 2024



Three Big Ones



A SPECIFIC *Salmonella Shigella* TABLE

Table 2A-2
Salmonella and *Shigella* spp.
CLSI M02 and CLSI M07

Table 2A-2. Zone Diameter and MIC Breakpoints for *Salmonella* and *Shigella* spp.

Testing Conditions		Routine QC Recommendations (see Tables 4A-1 and 5A-1 for acceptable QC ranges)
Medium:	Disk diffusion: MHA Broth dilution: CAMHB Agar dilution: MHA	<i>Escherichia coli</i> ATCC® 25922 <i>Pseudomonas aeruginosa</i> ATCC® 27853 (for carbapenems) <i>Staphylococcus aureus</i> ATCC® 25923 (for disk diffusion) or <i>S. aureus</i> ATCC® 29213 (for dilution methods) when testing azithromycin against <i>Salmonella enterica</i> ser. Typhi or <i>Shigella</i> spp.
Inoculum:	Broth culture method or colony suspension, equivalent to a 0.5 McFarland standard; positive blood culture broth for select antimicrobial agents with disk diffusion (see general comment [5])	
Incubation:	35°C ± 2°C; ambient air Disk diffusion: 16–18 hours Dilution methods: 16–20 hours	When a commercial test system is used for susceptibility testing, refer to the manufacturer's instructions for QC test recommendations and QC ranges.

- Ampicillin, fluoroquinolone, T/S for fecal; add 3° cephem for extra-intestinal isolates
- Testing is indicated for all *Shigella* spp. isolates

TABLE 2A-2

- Testable agents same as *Enterobacterales*

ampicillin

cefotaxime, ceftriaxone

ertapenem, imipenem, meropenem

tetracycline, doxycycline, minocycline

trimethoprim-sulfamethoxazole

chloramphenicol

- Azithromycin (different than *Enterobacterales*)

Only for *Shigella* spp., *Salmonella* serotype Typhi

Ciprofloxacin

TABLE 2A-2

Isolate	Broth Microdilution			Disk Diffusion		
	S	I	R	S	I	R
<i>Enterobacterales</i>	≤ 0.25	0.5	≥ 1	≥ 26	22-25	≤ 21
<i>Shigella</i> spp.	≤ 0.25	0.5	≥ 1	≥ 26	22-25	≤ 21
<i>Salmonella</i> spp.	≤ 0.06	0.12-0.5	≥ 1	≥ 31	21-30	≤ 20

Levofloxacin

Isolate	Broth Microdilution			Disk Diffusion		
	S	I	R	S	I	R
<i>Enterobacterales</i>	≤ 0.5	1	≥ 2	≥ 21	17-20	≤ 16
<i>Shigella</i> spp.	≤ 0.5	1	≥ 2	≥ 21	17-20	≤ 16
<i>Salmonella</i> spp.	≤ 0.12	0.25-1	≥ 2			

Ofloxacin/*Shigella* mirror *Enterobacterales*; no disk diffusion for *Salmonella*

BUH-BYE *Burkholderia cepacia* complex

Antimicrobial Agent	Disk Content	Interpretive Categories and Zone Diameter Breakpoints, nearest whole mm			Interpretive Categories and MIC Breakpoints, µg/mL			Comments
		S	I	R	S	I	R	
β-LACTAM COMBINATION AGENTS								
Ticarcillin-clavulanate*	–	–	–	–	≤ 16/2	32/2–64/2	≥ 128/2	
CEPHEMS (PARENTERAL) (Including cephalosporins I, II, III, and IV. Please refer to Glossary I.)								
Ceftazidime					≤ 8	16	≥ 32	
CARBAPENEMS								
Meropenem					≤ 4	8	≥ 16	
TETRACYCLINES								
Minocycline					≤ 4	8	≥ 16	
FLUOROQUINOLONES								
Levofloxacin	–	–	–	–	≤ 2	4	≥ 8	
FOLATE PATHWAY ANTAGONISTS								
Trimethoprim-sulfamethoxazole					≤ 2/38	–	≥ 4/76	
PHENICOLS								
Chloramphenicol*	–	–	–	–	≤ 8	16	≥ 32	(4) Not routinely reported on organisms isolated from the urinary tract.

CARBAPENEMS

- Resistance via:

Carbapenemases (direct hydrolysis of agent)
ESBL or AmpC + cell wall permeability defect

- New Delhi metallo- β -lactamase (NDM)

Hydrolyzes almost all traditional β -lactams

Not inhibited by ceftazidime-avibactam
imipenem-relebactam
meropenem-vaborbactam

Inhibited by aztreonam




FIG 5 Global distribution of metallo- β -lactamase-positive *Enterobacteriaceae* and *P. aeruginosa*, including NDM-type enzymes collected from 2012 to 2014 from surveillance. (Republished from reference 287).

THIS GETS COMPLICATED

- NDM isolates frequently harbor other β -lactamases

Able to hydrolyze aztreonam

Inhibited by avibactam

- Aztreonam and ceftazidime-avibactam (ATM-CZA)

Enhanced *in vitro* activity (next two slides)

Clinical efficacy against multi-drug- and

resistant to three or more classes

extensively drug-resistant

resistant to all but one or two classes

Enterobacterales (following two slides)

J Clin Microbiol. 61:e0164722; 2023

IN VITRO ASSESSMENTS

All Enterobacterales^a (N = 18 713)

Antimicrobials	MIC (mg/L)			%S CLSI	%S EUCAST ^b
	MIC ₅₀	MIC ₉₀	MIC range		
Aztreonam/avibactam ^c	0.03	0.25	0.015–128		99.9
Amikacin	2	8	0.25–128	95.2	92.8
Aztreonam	0.12	128	0.015–256	69.8	69.8
Cefepime	0.12	64	0.12–64	71.8	74.3
Ceftazidime	0.25	128	0.015–256	70.7	70.7
Colistin ^{d,e}	0.5	16	0.06–16	NA	97.1
Gentamicin	0.5	32	0.12–32	80.8	79.8
Imipenem	0.25	2	0.06–16	86.7	93.6
Levofloxacin	0.25	16	0.25–16	66.8	71.9
Meropenem	0.06	0.25	0.06–32	92.6	94.5
Piperacillin/tazobactam	4	128	0.12–128	81.3	75.5
Tigecycline	0.25	1	0.015–16	98.4 ^f	97.8 ^g

provisional MIC of $\leq 8 \mu\text{g/mL}$ used for susceptibility, based on pharmacokinetic/pharmacodynamic modeling

Organism	Breakpoint ($\mu\text{g/mL}$) for:				
	Aztreonam			Ceftazidime-avibactam	
	S	I	R	S	R
<i>Enterobacterales</i>	≤ 4	8	≥ 16	$\leq 8/4$	$\geq 16/4$
<i>Pseudomonas aeruginosa</i>	≤ 8	16	≥ 32	$\leq 8/4$	$\geq 16/4$
<i>Stenotrophomonas maltophilia</i>					

IN VITRO ASSESSMENTS

Table 5

Activity of aztreonam and aztreonam/avibactam (MIC in mg/L) against different enzyme variants and combinations for all Enterobacterales, 2019.

All Enterobacterales ^a (N = 18 713)	Drug	n	MIC (mg/L)			%S CLSI	%S EUCAST ^b
			MIC Range	MIC ₅₀	MIC ₉₀		
MBL positive ^c	Aztreonam	462	0.015–256	128	256	14.7	12.6
	Aztreonam/avibactam ^d		0.015–16	0.12	0.5		99.6
IMP ^e	Aztreonam	6	0.25–128	64	128	33.3	33.3
	Aztreonam/avibactam ^d		0.03–2	0.25	2		100.0
VIM ^f	Aztreonam	49	0.06–256	64	128	18.4	18.4
	Aztreonam/avibactam ^d		0.015–2	0.12	0.5		100.0
NDM ^g	Aztreonam	408	0.015–256	128	256	14.2	14.2
	Aztreonam/avibactam ^d		0.015–16	0.12	0.5		99.5
NDM-1	Aztreonam	270	0.015–256	128	256	14.4	14.4
	Aztreonam/avibactam ^d		0.015–4	0.12	0.5		100.0
NDM-5	Aztreonam	113	0.015–256	128	256	13.3	13.3
	Aztreonam/avibactam ^d		0.015–16	0.25	4		98.2
NDM-7	Aztreonam	17	0.03–256	128	256	23.5	23.5
	Aztreonam/avibactam ^d		0.03–0.5	0.12	0.5		100.0
IMP+VIM	Aztreonam	55	0.06–256	64	128	20	20
	Aztreonam/avibactam ^d		0.015–2	0.12	0.5		100.0
IMP+NDM	Aztreonam	414	0.015–256	128	256	14.5	14.5
	Aztreonam/avibactam ^d		0.015–16	0.12	0.5		99.5
NDM+VIM	Aztreonam	456	0.015–256	128	256	14.5	14.5
	Aztreonam/avibactam ^d		0.015–16	0.12	0.5		99.6
KPC positive ^h	Aztreonam	368	2–256	256	256	2.5	2.5
	Aztreonam/avibactam ^d		0.015–4	0.25	0.5		100.0
OXA positive ⁱ	Aztreonam	461	0.06–256	128	256	9.3	9.3
	Aztreonam/avibactam ^d		0.015–16	0.25	0.5		99.8
KPC+MBL positive	Aztreonam	820	0.015–256	128	256	9.4	9.4
	Aztreonam/avibactam ^d		0.015–16	0.25	0.5		99.8
OXA+MBL positive	Aztreonam	843	0.015–256	128	256	12.3	12.3
	Aztreonam/avibactam ^d		0.015–16	0.25	0.5		99.6
KPC+OXA+MBL positive	Aztreonam	1197	0.015–256	128	256	9.4	9.4
	Aztreonam/avibactam ^d		0.015–16	0.25	0.5		99.8

Efficacy of Ceftazidime-avibactam Plus Aztreonam in Patients With Bloodstream Infections Caused by Metallo- β -lactamase-Producing Enterobacterales

Marco Falcone,¹ George L. Daikos,² Giusy Tiseo,¹ Dimitrios Bassoulis,² Cesira Giordano,³ Valentina Galfo,¹ Alessandro Leonildi,³ Enrico Tagliaferri,¹ Simona Barnini,³ Spartaco Sani,⁴ Alessio Farcomeni,⁵ Lorenzo Ghiadoni,⁶ and Francesco Menichetti¹


¹Department of Clinical and Experimental Medicine, Infectious Diseases Unit, University of Pisa, Pisa, Italy, ²First Department of Medicine, School of Medicine, National and Kapodistrian University of Athens, Athens, Greece, ³Microbiology Unit, Azienda Ospedaliera Universitaria Pisana, Pisa, Italy, ⁴Infectious Disease Unit, Livorno Hospital, Livorno, Italy, ⁵Department of Economics and Finance, University of Rome "Tor Vergata," Rome, Italy, and ⁶Emergency Medicine Department, Azienda Ospedaliera Universitaria Pisana, University of Pisa, Pisa, Italy

- 102 bloodstream infections
 - 82 NDM; 20 VIM (carbapenemase)
 - 93 *Klebsiella pneumoniae*, 5 *Enterobacter* spp.
- 52 received ATM-CZA
- 50 received other active antibiotics (OAA)
- 27 with colistin

CLINICAL EFFICACY

Table 2. Targeted Antibiotic Regimens Administered in 102 Bloodstream Infections Due to Metallo- β -Lactamase-Producing Enterobacterales

Antibiotic Regimen	No. (%) (N = 102)	Mortality, No. (%)
CAZ-AVI + ATM ^a	52 (51)	10/52 (19.2)
OAA s		
Colistin-containing regimens	27 (26.5)	16/27 (59.3)
Colistin + fosfomycin + tigecycline	7	6/7
Colistin + fosfomycin	7	5/7
Colistin + meropenem	5	3/5
Colistin + ATM \pm piperacillin-tazobactam	4	1/4
Colistin + gentamicin	1	0/1
Colistin + cotrimoxazole	1	0/1
Colistin alone	2	1/2
Regimens not containing colistin	23 (22.5)	6/23 (26.1)
Tigecycline + aminoglycosides	8	2/8
Fosfomycin + aminoglycosides	5	0/5
Tigecycline + fosfomycin	2	2/2
Tigecycline + meropenem	1	0/1
ATM + aminoglycosides	4	1/4
ATM + fosfomycin	1	0/1
ATM alone	2	1/2





Number at risk	0	5	10	15	20	25	30
CAZ-AVI plus ATM	52	51	50	47	45	45	42
OAA	50	40	36	31	30	29	28

↓ 30d mortality rate $P = 0.007$
 ↓ d14 clinical failure $P = 0.002$
 shorter length of stay $P = 0.007$



Multicenter Evaluation of an MIC-Based Aztreonam and Ceftazidime-Avibactam Broth Disk Elution Test

Harley Harris,^a Lili Tao,^b Emily B. Jacobs,^a Yehudit Bergman,^a Ayomikun Adebayo,^a Tsigedera Tekle,^a Shawna Lewis,^a Ashley Dahlquist,^c Taylor C. Abbey,^c Eric Wenzler,^c  Romney Humphries,^b  Patricia J. Simner^a

- Broth disk elution method
 - 30 µg ATM
 - 30/20 µg CZA in 5 mL MH broth
 - 6/6/4 µg/mL ATM-CZA (growth/no growth)
- ~150 clinical isolates
 - metallo-β-lactamase *Enterobacterales*
 - carbapenem-resistant *P. aeruginosa*
 - Stenotrophomonas maltophilia*
- 97.9% categorical agreement vs. BMD; 2.4% ME

Table 3D. Aztreonam Plus Ceftazidime-Avibactam Broth Disk Elution Method¹

Due to limited therapeutic options, there may be a clinical need to assess the *in vitro* activity of the combination of aztreonam and ceftazidime-avibactam to guide therapeutic management of multidrug-resistant gram-negative bacterial infections, especially those caused by MBL producers.

The aztreonam plus ceftazidime-avibactam broth disk elution method was established with limited disk and/or media manufacturers and is considered provisional until additional data are evaluated by CLSI and shown to meet CLSI M23² guidance.

NOTE 1: Manufacturer-related issues were observed with different combinations of antimicrobial disks and CAMHB when the aztreonam plus ceftazidime-avibactam broth disk elution method was performed. QC of the method must be performed with every new lot or shipment of reagents to ensure the accuracy of results.

NOTE 2: Information in boldface type is new or modified since the previous edition.


Test	Aztreonam Plus Ceftazidime-Avibactam Broth Disk Elution
Organism group	Enterobacterales and <i>Stenotrophomonas maltophilia</i>
When to perform this test	Testing multidrug-resistant isolates, especially MBL producers
Test method	Tube dilution using aztreonam and ceftazidime-avibactam disks as the antimicrobial source
Medium	CAMHB (5-mL tubes)
Antimicrobial concentration	30-µg aztreonam disks 30/20-µg ceftazidime-avibactam disks Final concentration: 6 µg/mL aztreonam, 6 µg/mL ceftazidime, 4 µg/mL avibactam
Inoculum	1. Using a loop or swab, pick 3–5 colonies from a fresh (18–24 hours) nonselective agar plate and transfer to sterile saline (4–5 mL). 2. Adjust turbidity to equivalent of a 0.5 McFarland turbidity standard.

TABLE 3D

- Four 5-mL Mueller Hinton broth tubes; add disks

Mock	Ceftazidime-avibactam (1)
Aztreonam (1)	Aztreonam (1) + ceftazidime-avibactam (1)
- Vortex; allow 30-60 minutes for elution
- 25 μ L of 0.5 McFarland turbidity equivalent to all tubes
- Vortex at slow speed; ensure disks at bottom
- Incubate 16-20 hours in 33-35°C ambient air

TABLE 3D




susceptible to all
antimicrobial agents
evaluated

Escherichia coli ATCC 25922

CLSI M100-Ed34, 2024

TABLE 3D



not susceptible to ATM;
susceptible to CZA and
ATM-CZA

Klebsiella pneumoniae ATCC BAA-1705

CLSI M100-Ed34, 2024

TABLE 3D




not susceptible to ATM
or CZA; susceptible to
ATM-CZA

Klebsiella pneumoniae ATCC BAA-2146

CLSI M100-Ed34, 2024

TABLE 3D



not susceptible to any antimicrobial agents evaluated

this control is necessary due to manufacturer differences in disks and Mueller Hinton broth

Escherichia coli AR Bank #0348

CLSI M100-Ed34, 2024



Table 1



TABLE 1 IMPORTANT Δs

● Additions

Sulbactam-durlobactam *Acinetobacter* spp. (1B-2; tier 3)

New *Neisseria meningitidis* Table 1I

Penicillin/ampicillin Table 1J: Tier 1 Gram-positives

Tier 4 Gram-negatives

● Revisions

Wordsmithed comment about *Enterococcus*/penicillin

● Deletions

Stenotrophomonas maltophilia/ceftazidime (1B-4);
only cefiderocol remains



Table 2




TABLE 2A-1 IMPORTANT ADDITIONS

- Carbapenemase

Enterobacterales harboring OXA-48-like enzymes may test susceptible to meropenem-vaborbactam but have less clinical response; if OXA-48-like determinant or enzyme detected, suppress meropenem-vaborbactam result or report resistant

- More carbapenemase

Change cefepime S or SDD interpretations to resistant in isolates demonstrating carbapenemase

THE NON-FERMENTER TABLE 2

● Additions

Sulbactam-durlobactam MIC and disk diffusion breakpoints for *Acinetobacter* spp. (2B-2)

Trimethoprim-sulfamethoxazole should not be used for *Stenotrophomonas maltophilia* monotherapy


Organism	Method	Minocycline Previous			Minocycline New		
		S	I	R	S	I	R
<i>S. maltophilia</i>	BMD	≤ 4	8	≥ 16	≤ 1	2	≥ 4
	DD	≥ 19	15-18	≤ 14	≥ 26	21-25	≤ 20

● Deletions

All *Burkholderia cepacia* disk diffusion (2B-3)

CLSI M100-Ed34, 2024

GRAM-POSITIVE Δs



- Tedizolid additions

Disk diffusion for *S. aureus* (Table 2C)

Disk diffusion for *S. pyogenes* and *S. agalactiae* only (Table 2H-1)

Disk diffusion for *S. anginosus* group only (Table 2H-2)

- Linezolid

MIC confirmation no longer needed for *S. aureus* resistant via disk diffusion

Organism	Method	Linezolid Previous			Linezolid New		
		S	I	R	S	I	R
<i>S. aureus</i>	BMD	≤ 4		≥ 8	≤ 4		≥ 8
	DD	≥ 21		≤ 20	≥ 26	23-25	≤ 22



Table 3



AS A REMINDER...

Table 3A
Tests for ESBLs

Table 3B
CarbaNP Test for Suspected Carbapenemase Production



Table 3C
Modified Carbapenem Inactivation Methods

Table 3D
Aztreonam Plus Ceftazidime-Avibactam Broth Disk Elution Method

now Table 3E

Table 3D
Tests for Colistin Resistance for
Enterobacterales and *Pseudomonas aeruginosa*

now Table 3F-1,2,3,4

Table 3E
Test for Performing Disk Diffusion Directly From
Positive Blood Culture Broth

now Table 3G

Table 3F
Test for β -Lactamase Production
in *Staphylococcus* spp.

Table 3H
Oxacillin Salt Agar Test for Methicillin (Oxacillin) Resistance in *Staphylococcus aureus*

now Table 3I

Table 3H
Vancomycin Agar Screen for *Staphylococcus aureus*
and *Enterococcus* spp.

now Table 3J

Table 3I
Test for Inducible Clindamycin Resistance in *Staphylococcus* spp.,
Streptococcus pneumoniae, and *Streptococcus* spp. β -Hemolytic Group

now Table 3K

Table 3J
Test for High-Level Mupirocin Resistance in
Staphylococcus aureus

now Table 3L

Table 3K
Test for High-Level Aminoglycoside Resistance in
Enterococcus spp.



Test interpretation change from: positive, negative, indeterminate
to: positive, negative, inconclusive



Table 4



DISK DIFFUSION QC REVISIONS

<i>Staphylococcus aureus</i> ATCC 25923	tedizolid linezolid
<i>Staphylococcus aureus</i> ATCC 43300 ^a	cefoxitin

^a Listed as a supplemental strain (acceptable cefoxitin zone \leq 21 mm);
S. aureus ATCC 25923 also listed as QC strain (acceptable zone 23-29 mm)




Table 5




SOME MIC QC ADDITIONS/REVISIONS

<i>E. coli</i> ATCC 25922	upleganan aztreonam imipenem-funobactam
<i>E. coli</i> NCTC 13846	colistin QC alternative ^a polymyxin B QC alternative
<i>E. coli</i> ATCC BAA-3170	colistin QC alternative ^a
<i>K. pneumoniae</i> ATCC BAA-1705	imipenem-funobactam
<i>K. pneumoniae</i> ATCC 700603	aztreonam imipenem-funobactam
<i>P. aeruginosa</i> ATCC 27853	upleganan, colistin, imipenem-funobactam
<i>S. aureus</i> ATCC 43300	cefoxitin ($\geq 8 \mu\text{g/mL}$) oxacillin ($\geq 4 \mu\text{g/mL}$)
<i>S. aureus</i> ATCC 29213	exebacase

^a Colistin QC range has been deleted for *E. coli* ATCC 25922



Appendices



REVISIONS OF NOTE

● Appendix B (intrinsic resistance)

B1. Enterobacterales (Continued)

Antimicrobial Agent →	Ampicillin	Amoxicillin-clavulanate	Ampicillin-sulbactam	Ticarcillin	Cephalosporins I: Cefazolin, Cephalothin	Cephamycins: Cefoxitin, Cefotetan	Cephalosporins II: Cefuroxime	Imipenem	Tetracyclines	Tigecycline	Nitrofurantoin	Polymyxin B Colistin	Aminoglycosides
Organism ↓													
<i>Proteus vulgaris</i>	R				R		R	^d	R	R	R	R	
<i>Providencia rettgeri</i>	R	R			R			^d	R	R	R	R	
<i>Providencia stuartii</i>	R	R			R			^d	R	R	R	R	^e
<i>Raoultella</i> spp. ^f	R			R									
<i>Salmonella</i> and <i>Shigella</i> spp.	There is no intrinsic resistance to β -lactams in these organisms; refer to WARNING below for reporting.												
<i>Serratia marcescens</i>	R	R	R		R	R	R				R	R	
<i>Yersinia enterocolitica</i>	R	R		R	R								

Serratia marcescens/tobramycin removed

CLSI M100-Ed34, 2024

REVISIONS OF NOTE

● Appendix H (cefiderocol)



Figure H2-B




Figure H2-C

Aides in MIC determination; no trailing, no haze;
First well with button ≤ 1 mm is MIC

CLSI M100-Ed34, 2024



Thank you for your attention.
Have a better 2024.