

Foundational Training for the Entry-level Clinical Microbiologist Part Three: CLSI M35-A2 Document for Same-day Identification of Bacteria and Yeast (II)



Erik Munson
Wisconsin Clinical Laboratory Network
Laboratory Technical Advisory Group

The presenter states no conflict of interest and has no financial relationship
to disclose relevant to the content of this presentation.

OBJECTIVES

1. Apply basic bench-level skills, along with spot and rapid identification reagents, to the abbreviated identification of commonly-encountered patient isolates;
2. Recognize the utility of the CLSI M35-A2 guideline as an adjunct to routine operations in the clinical microbiology laboratory; and,
3. Review essential strategies in the triage of potential biothreat bacterial agents.



Life Without MALDI



Erik Munson
Wisconsin Clinical Laboratory Network
Laboratory Technical Advisory Group

The presenter states no conflict of interest and has no financial relationship to disclose relevant to the content of this presentation.

WE'LL DO BETTER NEXT TIME

- Foundational training

Hemolysis
Media

Odor
Significance

Colony morphology
Growth distribution
Cellular morphology

- M35-A2 Round One

Important reviews

- M35-A2 Round Two

Yeast Anaerobes
HACEKs Biothreats



(Re)Introduction to CLSI M35-A2

IMPETUS

- We've been doing this for years (odor, spot testing, *et cetera*); CLSI sought to standardize the practice
- Cost savings, patient benefits

Study	Turnaround Time <i>P</i>	Length of Stay <i>P</i> ⁴	Mortality Index <i>P</i> ⁴	Average Cost <i>P</i> ⁴	Notes
Doern <i>et al.</i> ¹	< 0.0005	NS ⁵	< 0.02	0.01	Rapid MicroScan product
Barenfanger <i>et al.</i> ²	0.001 ⁴	0.006	0.45	0.04	VITEK product
Kerremans <i>et al.</i> ³	< 0.0001	ND ⁶	0.21	ND	VITEK products

¹J Clin Microbiol 32:1757-1762; 1994

²J Clin Microbiol 37:1415-1418; 1999

³J Antimicrob Chemother 61:428-435; 2008

WARNINGS

- “Isolates to be tested should match the criteria required for proper identification.”
- Isolates conforming to M35-A2 identify organism with >95% accuracy; report without qualification
- Sniffing/wafting can be dangerous

Once mold colony is ruled out, opening of plates from non-invasive sources (i.e., urine, sputum) is common & relatively safe



CLSI M35-A2; 2008

OVERALL UTILITY I

	n	Ampicillin/Sulbactam	Piperacillin/Tazobactam	Imipenem	Meropenem	Aztreonam	Cefazolin	Ceftriaxone	Cefepime	Gentamicin	Amikacin	Ciprofloxacin	Levofloxacin	Trimethoprim/ Sulfamethoxazole
<i>Escherichia coli</i> ←	1595	58	96		100	96	89	95	96	93	100	80	80	79
<i>Klebsiella pneumoniae</i>	323	85	96		99	93	90	94	97	98	99	96	96	93
<i>Klebsiella oxytoca</i>	71	49	80		100	97	49	83	92	96	100	85	86	92
<i>Enterobacter cloacae</i>	89		87		97	90		88	98	98	100	94	94	86
<i>Enterobacter aerogenes</i>	38		71		100	84		87	97	100	100	89	89	97
<i>Citrobacter freundii</i>	30		93		100	93		93	100	97	100	96	97	90
<i>Proteus mirabilis</i> ←	155	84	99		100	100	87	100	100	95	100	81	85	81
<i>Serratia marcescens</i>	34				100	100		100	100	100	100	91	91	100
														Tobra
<i>Pseudomonas aeruginosa</i> ←	186		89	82	86	67			83	91		87	80	98
Organisms with extended-spectrum beta-lactamases (ESBL): 98														
<i>Haemophilus influenzae</i> beta-lactamase rate: 35.1%* ←														

OVERALL UTILITY II

	n	Penicillin	Methicillin	Erythromycin	Clindamycin	Vancomycin	Levofloxacin	Nitrofurantoin	Trimethoprim/sulfa	Tetracycline
methicillin-resistant <i>S. aureus</i> ← (45% MRSA rate)	255	0	0	10	69	100	41		100	91
methicillin-susceptible <i>S. aureus</i> ←	314	11	100	69	86	100	89		99	97
coagulase-negative <i>Staphylococcus</i> ←	221	13	48	28	50	100	65		66	91
<i>Enterococcus</i> species* ←	247-321	71				89		92		
viridans group <i>Streptococcus</i> ** ←	166	67		56	88	100				
<i>Streptococcus agalactiae</i> ** ←	167	100		39	54	100	98			
Beta-hemolytic streptococci other than group B** ←	40	100		83	85	100	100			
<i>Streptococcus pneumoniae</i> ** ←	81	PCN:	70.3% susceptible			Levofloxacin:		100.0% susceptible		
			18.5% intermediate					0.0% intermediate		
			11.1% resistant					0.0% resistant		

STOCK UP

Bile solubility reagent

Rapid spot CAMP

H₂O₂ (aerobes)

H₂O₂ (anaerobes)

H₂O₂ (*N. gonorrhoeae*)

Germ tube

Rapid hippurate hydrolysis

Spot indole

Indoxyl acetate

MUG test

Spot oxidase

δ-aminolevulinic acid

Rapid trehalose

Urea/phenylalanine deaminase


PYR, LAP

Latex agglutination (*S. aureus*)



SPOT REAGENTS

Reagent	Notes on Protocol	Read Time	Important Notes
Catalase	Smear then drop; different formulations	< 20 seconds	non-platinum loops; no blood agar
Indole	Different formulations	< 20 seconds	not on medium with dyes, not on Mueller-Hinton; indole diffuses
Oxidase	Slow-growers with no growth on MacConkey	< 10 seconds	no wires; no differential media



M35-A2 Rapid Identification



Candida albicans

- Presumptive identification

Budding yeast

- Additional tests for definitive identification

“Feet” in less than 48 hours OR
Germ tube-positive

- Notes

Not easily separated from
C. dubliniensis



Candida glabrata

- Presumptive identification

Small yeast in smear with no hyphae

Better growth on chocolate agar than blood agar

- Additional tests for definitive identification

Rapid trehalose-positive at 42°C OR

Evaluation of Growth Characteristics on Blood Agar and Eosin Methylene Blue Agar for the Identification of *Candida (Torulopsis) glabrata*

Martha J. Bale, Carol Yang, and
Michael A. Pfaller

Diagn Microbiol Infect Dis 28:65-67; 1997

TABLE 1 Growth Characteristics on Primary Culture Plates

	<i>C. glabrata</i>	Other Species
EMB growth > BAP ^a	92	0
EMB = or < BAP	10	267
Total	102	267

^a Colonies twice the size, or no growth on BAP.

Cryptococcus neoformans/gattii complexes

- Presumptive identification

Spherical pleomorphic budding yeast with no hyphae


- Additional tests for definitive identification

Urease-positive


Phenol oxidase-positive (caffeic acid disk)

- Notes

Phenotypic differentiation of complexes only by
canavanine-glycine bromthymol blue medium (5 d)



Gram stain
1000X total magnification







GRAM-POSITIVE COCCI

Parameter	viridans group <i>Streptococcus</i>	<i>Aerococcus</i> <i>viridans</i>	<i>Enterococcus</i>	<i>Aerococcus</i> <i>urinae</i>	<i>Pediococcus</i>
catalase	negative	negative	negative	negative	negative
morphology	pairs, chains	clusters, tetrad	pairs, chains	clusters, tetrad	clusters, tetrad
hemolysis	α or γ	α	γ or α	α	α
PYR	negative	positive	positive	negative	negative
LAP	positive	negative	positive	positive	positive
fun fact(s)	bile solubility- negative		watch out for α VRE	Grandma's UTI	vancomycin- resistant







Escherichia coli

- Presumptive identification

Gram-negative bacilli via Gram
or growth on selective medium

Oxidase-negative
Indole-positive

- Additional tests for definitive identification



Hemolytic OR lactose-positive AND PYR-negative
OR lactose-negative AND MUG-positive


- Notes

Isolate must be growing as large colonies

Isolate cannot be derived from gastrointestinal source
(*Shigella* spp.)

Escherichia coli





Let's (finally) Finish this Up



Neisseria meningitidis

- Presumptive identification

Gram-negative diplococci

Oxidase-positive

Glistening non-hemolytic growth on blood agar

- Additional tests for definitive identification

γ -glutamyl-aminopeptidase-positive

- Notes

Highly infectious; work in biological safety cabinet

Neisseria lactamica

- Presumptive identification

Gram-negative diplococci

Oxidase-positive

Graying non-hemolytic growth on blood agar

- Additional tests for definitive identification

β -galactosidase-positive

- Notes

Neisseria gonorrhoeae

- Presumptive identification

Gram-negative diplococci Oxidase-positive
Strongly catalase-positive with 30% H₂O₂

- Additional tests for definitive identification

β-galactosidase-negative AND
γ-glutamyl-aminopeptidase-negative

- Notes

Growth on *N. gonorrhoeae*-selective medium OR
No growth on Mueller-Hinton or tryptic soy media

Moraxella catarrhalis

- Presumptive identification

Gram-negative diplococci

Hockey puck on blood agar

Oxidase-positive

Catalase-positive

- Additional tests for definitive identification

Butyrate esterase-positive OR

Indoxyl acetate-positive

- Notes

Broth microdilution and disk diffusion susceptibility testing guidelines described in M45-A2 (Mueller-Hinton media)

Moraxella catarrhalis

- Presumptive identification

Gram-negative diplococci
Hockey puck on blood agar

Oxidase-positive
Catalase-positive




- Additional tests for definitive identification

Butyrate esterase-positive
Indoxyl acetate-positive

LIMITATIONS/additional factors

Most other *Moraxella* spp. are positive for butyrate esterase, but are coccobacilli rather than diplococci

RUSH



SHOW DON'T TELL



GRAM-NEGATIVE COCCOBACILLI

Genus	Growth Rate	Blood agar	Chocolate agar	MacConkey agar	Fun Fact(s)
<i>Haemophilus</i>	overnight	no growth	growth	no growth	M35-A2
<i>Pasteurella</i>	overnight	growth	growth	no growth	sick TSI, penicillin
<i>Gardnerella</i>	overnight	growth	growth	no growth	teeny tiny
<i>Brucella</i> (nasty one)	two days	growth	growth	no growth	BSC, rule-out
<i>Francisella</i>	two days	no growth*	growth	no growth	BSC, rule-out
<i>Acinetobacter</i>	overnight	growth	growth	growth	non-motile, K/NC



GRAM-NEGATIVE COCCOBACILLI

Genus	Growth Rate	Blood agar	Chocolate agar	MacConkey agar	Fun Fact(s)
<i>Haemophilus</i>	overnight	no growth	growth	no growth	M35-A2
<i>Pasteurella</i>	overnight	growth	growth	no growth	sick TSI, penicillin
<i>Gardnerella</i>	overnight	growth	growth	no growth	teeny tiny
<i>Brucella</i> (nasty one)	two days	growth	growth	no growth	BSC, rule-out
<i>Francisella</i>	two days	no growth*	growth	no growth	BSC, rule-out
<i>Acinetobacter</i>	overnight	growth	growth	growth	non-motile, K/NC
HACEKs	overnight	WE	WILL	SEE	SOON



Cardiobacterium hominis

- Presumptive identification

Pleomorphic, thin GNR
No growth on MacConkey

Oxidase-positive
Catalase-negative

- Additional tests for definitive identification

Indole-positive
Non-hemolytic on blood agar; may pit

- Notes

Isolate must be derived from blood culture
Confirm with negative nitrate test
CLSI M45; CAMHB with LHB for all of the ACEKs

Cardiobacterium hominis

- Presumptive identification

Pleomorphic, thin GNR
No growth on MacConkey

Oxidase-positive
Catalase-negative

- Additional tests for definitive identification

Indole-positive
Non-hemolytic on blood agar; may pit

LIMITATIONS/additional factors

Pasteurella bettyae has same rapid biochemical profile as *C. hominis*, but is not known to cause endocarditis

Rosette-forming, thin GNR (with above biochemicals) in adult blood cultures

Eikenella corrodens

- Presumptive identification

Small GNR; pits blood, chocolate	Oxidase-positive
No growth on MacConkey	Catalase-negative

- Additional tests for definitive identification

Indole-negative	Distinct odor of bleach
Non-hemolytic on blood agar	

- Notes

Capnophilic

The only MacConkey-negative, catalase-negative, oxidase-positive GNR that is ornithine-positive

Kingella kingae

- Presumptive identification

Gram-negative coccoid bacillus
No growth on MacConkey

Oxidase-positive
Catalase-negative


- Additional tests for definitive identification

Hemolytic colonies on blood agar

- Notes

Joint fluids, blood cultures, other sterile body sites
Early growth confused for β -hemolytic streptococci


MISSOURI



ANAEROBIC GRAM-NEGATIVES

Table 4. Abbreviated Identification of Anaerobic Gram-Negative Bacteria

Identification	aBAP/LKV Colony Morphology	Cell Morphology	Bacteroides Bile Esculin Colony Morphology	Indole Reaction
<i>Bacteroides fragilis</i> group	Large, convex	Regular	Large, convex, gray-black	Not Done
<i>Bacteroides ureolyticus</i>	Translucent, pitting the agar (catalase-negative with 15% H ₂ O ₂)	Tiny rods, coccobacilli	No growth	Negative
<i>Bilophila wadsworthia</i>	Tiny, translucent (catalase +++ with 15% H ₂ O ₂)	Regular to filaments	Translucent with black center at 72 hours	Negative
<i>Fusobacterium nucleatum</i>	Opalescent, breadcrumb	Fusiform, thin pointed	No growth	Positive
<i>Porphyromonas</i> spp.	Small, translucent or opaque, fluoresce brick-red on aBAP	Tiny coccobacilli	No growth	Positive
<i>Prevotella intermedia</i>	Small, translucent or opaque, fluoresce brick-red on LKV or aBAP	Tiny coccobacilli	No growth	Positive
<i>Prevotella</i> spp.				Negative
<i>Veillonella</i> spp.	Small, transparent or opaque, fluoresce brick-red on aBAP	Tiny diplococci	No growth	Negative



ANAEROBIC GRAM-POSITIVES

Table 5. Abbreviated Identification of Anaerobic Gram-Positive Bacteria

Identification	aBAP Colony Morphology	Cell Morphology	Indole Reaction
<i>C. difficile</i>	Large, flat colonies; barnyard (cow manure) smell; chartreuse fluorescence	Thin rods, rare spores	Negative
<i>C. perfringens</i>	Large, irregular-shaped, double zone beta-hemolysis	Boxcar, large, square rods	Not done
<i>C. septicum</i>	Smoothly swarming	Thin rods, subterminal spores	Negative
<i>C. sordellii</i>	Very large, lobate, irregular, flat	Thin rods, subterminal spores	Positive
<i>C. tetani</i>	Smoothly swarming but slow growing	Swollen terminal spores	Positive
<i>Peptostreptococcus</i> spp.	Small, peaked, circular	Cocci, pairs and chains	Not done
<i>Propionibacterium acnes</i>	Small, opaque, enamel-white, circular (catalase+ with 15% H ₂ O ₂)	Coryneform rods	Positive





CLSI M35-A2, BIOTHRREAT AGENTS

- Changes in recommendations for handling cultures on the basis of unsuspected exposure to agents of bioterrorism
- Colonial growth examined in biological safety cabinet (while wearing gloves) until agents of bioterrorism or highly-pathogenic agents ruled out

Blood cultures

CSF cultures

Lymph node cultures

MORE...

- Unidentified Gram-negative, Gram-variable bacillus, or Gram-negative coccobacillus that only grows on blood and chocolate agar only (not MacConkey) is handled with extreme caution until rule-out
- Gram staining and wet mount preparation takes place in biological safety cabinet; wear gloves

THIS IS SERIOUS

- Procedures known to create aerosols (catalase, others) should be confined to biological safety cabinet with extra care or avoided all together until rule-out
- Automated systems may pose danger with respect to these organisms; may not generate accurate identification or susceptibility data in the first place

Should not be performed until bioterrorism agent possibility has been eliminated

J Hosp Infect 22:159-162; 1992

Diagn Microbiol Infect Dis 60:241-246; 2008

Brucella spp.

- Presumptive identification


Tiny Gram-negative coccobacillus Oxidase-positive
No growth on MacConkey Catalase-positive

- Additional tests for definitive identification (**R/O**)

Urease-positive Indole-negative
Non-hemolytic on blood agar

- Notes

Phenylalanine-negative
Sterile tissues and fluids
Growth on CNA agar



Francisella tularensis

- Presumptive identification


Tiny GNR or coccobacillus (sand) Oxidase-negative
Slow on chocolate; no satellite Catalase-neg/wk

- Additional tests for definitive identification (**R/O**)




β -lactamase-positive

- Notes

L-cysteine
On the rise?



STAY TUNED



Acknowledgments

Jim Hermanson

Kacy Fisher

Joshua Kropp

Alana Sterkel, Ph.D.

David Uselman

Nicholas Zumwalde, Ph.D.

Wisconsin Diagnostic Laboratories

Wisconsin State Laboratory of Hygiene

Wisconsin Clinical Laboratory Network

Laboratory Technical Advisory Group