







Systemic Mycoses

Dave Warshauer, PhD, D(ABMM) Deputy Director, Communicable Diseases

david.warshauer@slh.wisc.edu





Systemic Mycoses

- Histoplasma capsulatum
- Blastomyces dermatitidis
- Coccidioides immitis
- Coccidioides posadasii
- Paracoccidioides brasiliensis

- Penicillium marneffei
- Sporothrix schenckii
- Aspergillus species
- Emmonsia species





Dimorphism

- Majority systemic mycoses are caused by dimorphic fungi
 - Molds in their normal soil environment, and on laboratory media when kept at 25-30°C.
 - Yeast when the temperature is raised to 37°C (as in the human host).
 - Nutritional factors such as certain amino acids can also enter into dimorphism, but are generally less important than temperature.



Histoplasmosa capsulatum var. capsulatum



Nethealthbook.com





Histoplasmosis

- Occurs throughout the world
- Endemic areas
 - Mississippi and Ohio River Valleys in the US.
 - Mexico, Central and South America
- Once thought to be a highly lethal form of pneumonia with up to 90% mortality
- Now known to be a rather common infection in endemic areas.





Histoplasma habitat

- Soil saprobe
- Loves droppings from bats or birds.
 - Requires high levels of creatinine and nitrogen
- Birds not susceptible to infection with Histoplasma
 - · Likely related to their high body temperature
 - 40-42C (104-108F)
- Occupational risk for people working with chickens.
- Clearing Starling roosts has been associated with large outbreaks of fatal infections.



Histoplasma capsulatum Disease

- Usually self limiting flu-like illness and does not require medical intervention.
 - 90% asymptomatic
 - 4:1 male predominance for clinical disease
 - Infants and young children more likely symptomatic
 - Chronic pulmonary disease
 - Disseminated disease in immunocompromised



Pathobiology

In human infection---small yeast 2-5 µm in diameter.

- predominantly in macrophages.

- Nonactivated macrophages do not effectively kill <u>H. capsulatum</u> and can actually spread the disease.
 - Can multiply intracellularly, kill the phagocyte, and infect additional cells



Pathobiology

- Granulomatous lesions in the lungs
 - Very similar to tuberculosis lesions
 - Lymphocytes, macrophages, Langhans' giant cells
- Severity of infection directly proportional to the number of conidiospores inhaled.
 - Miliary lesions when large numbers of the spores are inhaled.
- In most infections the cure is spontaneous and lasting immunity occurs.
- Histoplasma may remain viable and recurrence possible with decrease in CMI



Histoplasmosis

Diffuse pneumonic // histoplasmosis

•Radiologic variations

Calcified miliary histoplasmosis



Diffuse pneumonic lesions throughout both lungs, representing acute or epidemic histoplasmosis

Millary histoplasmosis





Many small parenchymal and hilar areas of calcification in both lungs; classic appearance of heated histoplasmosis Bilateral infiltrates with cavitation in 1 upper lobe. This chronic progressive, cavitary form of histoplasmosis appears identical to tuberculosis operation





Calcification







Immunity and Treatment

- Immunity dependent on CMI.
 - Antibody is of little importance
- Healing of lesions leads to calcified granulomas similar to that seen in tuberculosis.
 - Old calcified nodules on chest x-ray not uncommon
- Treatment reserved for life-threatening infections
 - Amphotericin B
 - Itraconazole





Laboratory Diagnosis

- Histology
 - GMS stain
 - Wright stain of blood or bone marrow
- Direct Microscopic Exam
 - KOH, Calcofluor---2-5um yeast
- Culture
 - <u>Enriched media</u> (BHI with Blood, Yeast Extract Phosphate, Inhibitory Mold Agar)
 - 2-4 weeks at 30C
- Antigen Detection-----Urine EIA (miravistalabs.com)
- Serology
 - Complement fixation
 - EIA
 - Immunodiffusion
- Real-time PCR



Specimens for Fungal Infections

- Respiratory secretions, tissues, blood, CSF, other body fluids
- Discourage swabs
- Transport at room temp
 - Specimens with endogenous flora, refrigerate if >2hr delay
- Blood--Lysis centrifugation, BACTEC MYCO/F Lytic, or BacT ALERT MB
- CSF—large volume (10-20ml)
 - Centrifuge 2000g, 10 minutes
 - Inoculate pellet





Specimens for Fungal Infections

- Urine and other body fluids
 - Centrifuge 2000g, 10 minutes
 - Plate pellet
- Mince tissue, do not grind
 - Place 3-4 pieces on plate and press into agar
 - Exception for Histoplasma—want to grind to release intracellular organisms



Media for Primary Isolation of Systemic Fungi

- Non-inhibitory media
 - Sabouraud's dextrose agar
 - Potato Flake Agar
 - Potato Dextrose Agar
- Selective Media
 - Mycobiotic or Mycosel agar--cyclohexamide and chloramphenicol



Media for Primary Isolation of Systemic Fungi (2)

- Enriched media w/ or w/o antibiotics
 - Inhibitory mold agar---chloramphenicol and cyclohexamide
 - BHI with sheep blood w/wo antibiotics
 - Yeast extract phosphate agar with ammonia
- Incubate plates or tubes at 30°C or 25°C
 - Hold 4 weeks
- For Blood Cultures
 - Lysis Centrifugation
 - BACTEC MYCO/F or BacT ALERT MB



Silver Stain





















Culture



Histo at 3 weeks, 30C





Histoplasma capsulatum



totallyfreeimages.com



totallyfreeimages.com





H. capsulatum conidia 30C Incubation



WELLE WISCONSIN STATE LABORATORY OF HYGIENE



H. capsulatum







H. capsulatum—Yeast phase







Differentiation from other Fungi

- Must differentiate from *Sepedonium* and *Chrysosporium* species that produce tuberculate macroconidia
 - More rapid growing
 - Not dimorphic
 - Usually will not grow in the presence of cycloheximide
 - Distinguish using DNA probe





Nucleic Acid Probe Identification

- GenProbe® Assay
 - Rapid
 - Chemiluminescent assay using labeled probes specific for each agent
 - Labeled DNA probe hybridizes with rRNA of the fungus
 - Available for <u>H. capsulatum</u>, <u>Blastomyces</u>
 <u>dermatitidis</u>, and <u>C. immitis</u>













Blastomycosis

- Blastomyces dermatitidis
- Agent of North American
 Blastomycosis,
- Geographical distribution is similar to <u>H. capsulatum</u>
- More common in Wisconsin than <u>H.</u> <u>capsulatum</u>.





Kurt Reed et. al PLOSone 3(4): e2034, 2008



WELLE WISCONSIN STATE LABORATORY OF HYGIENE

A

В



Blastomycosis

- The epidemiology is poorly understood
 - Lack of a good skin test reagent
 - Ecologic niche not well established
- Difficult to recover from the soil in endemic areas.
- Eagle River, Wisconsin outbreak 1985
 - First time <u>Blastomyces</u> isolated from the environment at the site of an outbreak
 - Isolated from soil containing decayed vegetative matter and from decomposed wood.



Clinical Manifestations

- Two clinical presentations
 - A primary cutaneous infection which usually remains localized to one area of the body
 - May indicate systemic disease
 - Primary pulmonary infection with possible secondary dissemination.
 - 30-45 day incubation
 - Mimics flu progressing to cough, weight loss, chest pain, low grade fever
 - 75% with isolated pulmonary disease
 - Infection may involve any organ
 - Secondary cutaneous infection
 - Asymptomatic in >50% of those infected



Systemic Disease

- Common sites of infection in systemic disease
 - Bones---long bones, ribs, vertebrae
 - Joints
 - Genitourinary tract---prostate, epididymis
 - CNS-----common in AIDS (40%), uncommon in immunocompetent (<5%)



Cutaneous Form

- A chronic suppurative granulomatous lesion.
- The presence of epithelial microabscesses and characteristic yeasts in the tissues is considered diagnostic.
- It is important to obtain urine and sputum samples from a patient with cutaneous blastomycosis since systemic spread may occur.





Histology

Body fluids or tissue specimens

- Look for the characteristic yeast form.

- Large (8-15 µm) and thick walled.
- The wall is prominent; "doubly refractile" on bright field microscopy.
- A single daughter cell (bud) is present with a broad connection between the two cells (BROAD-BASED BUDDING).



















Direct Exam



KOH Prep

Calcofluor White





KOH Exam









Calcofluor White







Gram Stain







Culture Characteristics

- On culture:
 - Slow growing gray/white mold
 - Delicate, septate hyphae
 - Conidia usually absent on bloodcontaining media. May be sparse on PDA and SAB
 - "Lollypop" conidiation





Blastomyces dermatitidis













Blastomyces Mold Phase

- This form of conidia is also found in such fungi as Chrysosporium sp., Pseudallescheria boydii (Scedosporium), and various Trichophyton sp.
- Differentiation from these other species can be made by the following characteristics:
 - Slower growth
 - Growth in the presence of cycloheximide
 - Dimorphism
 - Nucleic acid probes





Blastomyces Yeast Phase







Blastomyces Conversion







WETER WISCONSIN STATE LABORATORY OF HYGIENE



Coccidioides immitis and C. posadasii







HEAD OF DOMINGO EZCURRA FIRST CASE OF COCCIDIOIDOMYCOSIS





PREVALENCE OF COCCIDIODIN SENSITIVITY IN YOUNG ADULTS







Coccidiodes immitis

- Coccidiomycosis is sometimes known as "<u>San Joaquin</u> <u>Valley fever</u>". Up to 95% of the residents of the endemic area are skin test positive (coccidioidin test positive)
- Lower Sonoran Life Zone
 - Arid climate, hot summers, few winter freezes, low altitude, alkaline soil, sparse flora
 - Drought followed by heavy rains---Increased infections
 - 100,000 infected annually in U.S.
- Variety of animals infected
 - Positive cultures around rodent burrows
- Archaeology students discover new "infected" sites



Clinical Manifestations

- The primary disease is pulmonary, secondary to inhalation of small numbers of arthrospores
- Usually resolves spontaneously as an influenza-like infection.
 - 60% asymptomatic
 - 40% influenza-like illness, LRI or systemic illness
 - Cough, sputum, chest pain, malaise, fever, chills, night sweats, arthralgias, anorexia



Clinical Manifestations

- In a minority of cases a more chronic pulmonary infection occurs
 - Granulomatous lesions of the lung
 - Can lead to cavitation
- In rare cases (0.5%) dissemination occurs which can lead to rapidly fatal results.
- Reactivation infection occurs



Coccidiomycosis

Lesions variable: Papules Pustules Plaques Nodules Ulcers Abscesses









Histology

- Histological examination useful in confirming a diagnosis.
 - Spherules
 - 10-60 μm in diameter, but they may be as big as 200 $\mu m.$
- Immature spherules can be similar in size to the large yeast cells of <u>*B. dermatitidis*</u>
- The spherules contain <u>endospores</u> 2-5 µm in diameter
 - Similar in size to Histoplasma capsulatum
 - Will not see budding









Culture

- Culture:
 - The organism grows fairly rapidly.
 - Visible growth on Sabouraud's agar within a few days.





Microscopic Features

• Arthrospores:

- Formed by fragmentation of hyphae
- Very thick walled. Provides them with resistance to drying.
- "Barrel-Shaped"
 - As a culture ages on media the entire hyphal mass may fragment and form arthrospores.
- <u>Spores are extremely infectious Handle with</u> <u>extreme care</u>.
- **Remember**: Arthrospores can be made by other fungi.
 - Malbranchea sp., Gymnoascus uncinatus, Auxarthron sp.
 - <u>Geotrichum</u> and <u>Trichosporon</u> can also form arthrospores
- Confirm the identification nucleic acid probe test



Coccidioides---Arthrospores





Select Agent Regulations

- Report to CDC within 7 days of ID
 - Responsibility of lab confirming ID
 - Select Agent APHIS/CDC Form 4
- Secure against loss, theft, or release
- Destroy all subcultures and specimens
- Good News
 - Proposed to remove Coccidioides from SA list















WEILE WISCONSIN STATE LABORATORY OF HYGIENE



