

The Changing Landscape of Stool Parasite Diagnosis and Surveillance

Ann Valley- Advanced Microbiologist
Tim Monson- Microbiologist Supervisor
Communicable Disease Division
Wisconsin State Laboratory of Hygiene

WCLN Audioconference- Feb 11, 2015

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN


Objectives

- Discuss intestinal parasites seen in WI and some resources available to aid in their diagnosis
- Discuss continued use of traditional diagnostic methods for parasite detection
- Describe multi-target diagnostic assays and rapid antigen tests that are currently available for parasite detection and their effect on reporting and disease surveillance

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

Let's talk about coccidian parasites:


- *Cryptosporidia*
- *Cyclospora cayentanensis*
- *Cystoisospora*



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

Cryptosporidium spp.

- Most common cause of waterborne disease in the US.
- Most common parasite infection in WI
- Spread through contaminated food or water or direct or indirect contact with human or animal feces.
- Produces watery diarrhea 2-10 days post exposure which may last up to 2 weeks



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN


Cryptosporidium spp.





WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

Cryptosporidium spp. Community Pool Outbreak 2013



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

7

Cryptosporidium spp.

- Each person has an average 0.14 grams of fecal material on their perianal surface if they do not take a pre-swim shower with soap.
- A single diarrheal accident can introduce 10^7 - 10^8 *Cryptosporidium* oocysts into the water- enough to cause infection with a single mouthful of pool water.
- Crypto survives up to 11 days at chlorine concentrations found in most pools.

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

8

Cyclospora cayetanensis

- Found in tropical and subtropical regions
- Outbreaks in the US are associated with contaminated fresh produce.
- Incubation: ~ 1 week
- Watery diarrhea, cramping, low-grade fever
- Symptoms may last for several weeks
- Treatment with SXT



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

9

Cyclospora cayetanensis

Summer 2013

- June 15-29: 86 Nebraska and 153 Iowa cases associated with salad consumption at 2 national restaurant chains
- July-August: 278 Texas cases associated with cilantro consumption at local restaurant.

<http://www.fda.gov/food/recalls/outbreakemergencies/outbreaks/ucm361637.htm>



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

10


Cystoisospora belli

- Worldwide distribution
- More common in immunocompromised individuals
- Institutional outbreaks have occurred in the US
- Presents with non-bloody diarrhea which can last for weeks; more severe in children and elderly

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN


11


The Flagellates: *Giardia* and *D. fragilis*




WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

12









Giardia spp.

- Global distribution
- Diarrhea, malabsorption (1-2 weeks)
- Risk factors:
 - Travel to endemic areas
 - Backpackers, campers
 - Consumption of contaminated drinking water or water from lakes, rivers
 - child-care workers




WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN
13




Giardia spp.

Giardiasis associated with transient well contamination:

August 2014




WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN
14




Dientamoeba fragilis

- Worldwide distribution
- Diarrhea and abdominal pain.
- Possible association with pinworm infections
- Fragile and may not survive well outside host
- Increased risk for those living under poor sanitary conditions





WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN
15




E.histo/dispar and *B. hominis*

- *Entamoeba histolytica/dispar*
 - More common in tropical areas
 - Illness can range from asymptomatic to amebic dysentery.
 - Complications may include organ invasion
 - Important to differentiate species


- *Blastocystis hominis*
 - Clinical significance is questionable.
 - Disease may be associated with parasite load or with specific subtypes.
 - Suggestion of correlation of *B. hominis* infection with Irritable Bowel Disease (IBS)
 - Traveler's diarrhea.

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN
16




Balantidium coli




- Infection occurs when the cysts are ingested by eating contaminated food or water.
- Approximately 1% of people worldwide are infected.
- Occurs mainly in developing countries.
- Higher risk among pig farmers
- Diarrhea, weight loss, dysentery

Microsporidium






- Worldwide distribution
- 1200 species (15 known to infect humans)
- More common in immunocompromised patients (HIV)
- Many clinical manifestations including diarrhea, corneal and muscular infection

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN
17



WORMS!!!!

- Nematodes
- Cestodes
- Trematodes






WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN
18

Nematodes



Whipworm



Ascaris



Hookworm


WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

19

Nematodes

Enterobius vermicularis: Pinworm


- Most common worm infection in the US
- Occurs in young children and may spread to family members
- Eggs may survive 2-3 weeks on clothing and surfaces.



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

20

Diphyllobothrium latum (Freshwater fish tapeworm)



- Largest tapeworm that can infect humans
- Consumption of raw or undercooked fish
- Predominantly in Northern hemisphere
- Most cases asymptomatic but may cause diarrhea, intestinal obstruction, gall bladder disease

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

21

Taenia spp.

- Consumption of undercooked beef and pork
- Digestive issues with weight loss
- Cysticercosis (*T. solium*)
- Recognized as one of the 5 Neglected Parasitic Diseases in the US by CDC.



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

22

Liver Flukes




- Ingestion of the of metacercarial stage of the parasite.
- Adult flukes develop in ~3 months
- Intestinal obstruction and accumulation of abdominal fluids due to the blockage of bile ducts.


WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

23

Trichrome Stain





E. histolytica/dispar



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

24

Trichrome Stain







E. histolytica/dispar

Entamoeba coli

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN 25

Modified Acid Fast and Hot Safranin

Cryptosporidium

Cyclospora

Cystoisospora belli

Images from CDC DPDx




WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN 26

Wet Mount

Protozoa						
<i>Cystidia polydiplopori</i>	<i>Entamoeba histolytica</i>	<i>Trichomonas vaginalis</i>	<i>Ascaris lumbricoides</i>	<i>Ascaris lumbricoides</i>	<i>Hookworms</i>	<i>Trichostrongylus axei</i>
Cestodes						
<i>Taenia solium</i>	<i>Hymenolepis nana</i>	<i>Hymenolepis diminuta</i>	<i>Dipylidobothrium latum</i>	<i>Dipylidobothrium latum</i>		
Scale 100 µm						

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN 27



Wet Mount

- *Giardia* trophozoite stained using Trichrome
- *Giardia* cyst stained using Trichrome
- *Giardia* cyst from wet mount using Iodine

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN 28

Epifluorescence



Cyclospora cayetanensis


Cystoisospora belli

Images from CDC DPDx

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN 29

Identification of worms



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN 30

Non-Traditional Parasite Dx Tests

- Serology
- Microplate EIA
- Direct Fluorescent Antibody
- Rapid Cartridge Assays
 - Optical immunoassay
 - Lateral flow immunoassay
- Polymerase Chain Reaction (PCR)- Single target
- Multi-target Assays
 - PCR- based
 - Film array
 - Bead-based technology

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

31

Serology- Stool Parasites

- No serology tests available at WSLH for diagnosis of stool parasite infections
- CDC offers a limited menu of test options for stool parasite detection in serum
 - Antibody detection
- Reference laboratories might offer serology testing for the more common (to the U.S.) systemic stool parasites

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

32

Serology Testing- CDC Antibody Detection

Parasite	Test	Specimen
<i>E. histolytica</i>	EIA	Serum
<i>Baylisascaris</i>	Immunoblot	Serum or CSF
Cysticercosis (<i>T. solium</i>)	Immunoblot	Serum or CSF
<i>Echinococcus</i>	EIA/ Blot	Serum
<i>Paragonimus</i>	Blot	Serum
<i>Schistosoma</i>	Fast-ELISA/ Blot	Serum
<i>Strongyloides</i>	EIA	Serum
<i>Toxocara</i>	EIA	Serum or Vitreous
<i>Trichinella</i>	EIA	Serum

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

33

Antigen Detection

- Microplate EIA
 - High throughput
 - May cost more in long run if kits not used up
 - Generally very good sensitivity and specificity
- DFA
 - Very good sensitivity and specificity
 - Requires fluorescent microscope
- Rapid cartridge assays
 - Easy to perform
 - Assess performance, PPV/NPV before implementation

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

34

<http://www.cdc.gov/dpdx/diagnosticProcedures/stool/antigendetection.html>

Organism	Kit name	Manufacturer - distributor*	Type of Test*
Cryptosporidium spp.	Crypto CELISA	Cellabs	EIA
	PARA-TECT™ Cryptosporidium Antigen 90	Medical Chemical Corporation	EIA
	ProSpecT Rapid	Biorad	EIA
	ProSpecT	Biorad	EIA
	Cryptosporidium	TrichLab	EIA
	Cryptosporidium	Wampole	EIA
	Crypto GEL	Cellabs	EIA
Cryptosporidium spp./Giardia duodenalis	SPect Crypto	Biorad	Rapid
	PARA-TECT™ Cryptosporidium/Giardia DFA 95	Medical Chemical Corporation	DFA
	Meridian	Meridian	DFA
	ProSpecT	Biorad	EIA
	Crypto-Giardia CEL	Cellabs	EIA
	Colo-PAC	Becton Dickinson	Rapid
	ImmunoCard STAT®	Meridian	Rapid
Cryptosporidium spp./Giardia duodenalis/Entamoeba histolytica/Shigella	SPect	Biorad	Rapid
	Trigeo	BioMérieux	Rapid
	Entamoeba CELISA	Cellabs	EIA
	E. histolytica	Wampole	EIA
	E. histolytica II	TrichLab	EIA
	ProSpecT	Biorad	EIA
	Giardia CELISA	Cellabs	EIA
Giardia duodenalis	PARA-TECT™ Giardia Antigen 90	Medical Chemical Corporation	EIA
	ProSpecT	Biorad	EIA
	Giardia II	TrichLab	EIA
	Giardia	Wampole	EIA
	Giardia EIA	Antibodies, Inc.	EIA
	Giardia CEL	Cellabs	EIA
	ProSpecT	Biorad	Rapid
Wuchereria bancrofti	Simple-Screen Giardia	Medical Chemical Corporation	Rapid
	Fluoracut CELISA	Cellabs	EIA

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

35

Microplate EIA

- Many commercial kits on market for *Crypto* and *Giardia*
- Few kits available for *E. histolytica*/ *dispar* (only Techlab EIA assay will differentiate *histolytica*/ *dispar*)



Google Images

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

36

Direct Fluorescent Antibody

- Many consider this a gold standard test
- Stool concentration recommended prior to testing to increase sensitivity*
- Read using different wavelength filter than that used for epifluorescence
- Easily adapted to a lower throughput lab
- Generally easy to read and interpret
 - Background minimal with brightly fluorescing cysts or oocysts
 - Note size and morphology of cysts/ oocysts

* Most antigen detection tests do not recommend concentration prior to testing


WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

37

Direct Fluorescent Antibody


Cryptosporidium oocysts

image from CDC



Giardia cysts

Image from CDC



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

38


Rapid Cartridge Assays

- Concentration not recommended prior to use
- Very easy to use; Must strictly adhere to the procedure times in the package insert
- In general, stool antigen flows across a membrane containing antibody against the targeted parasite(s); Ab-Ag reaction leads to an immunochromatographic (colorimetric) reaction within a defined time period
- Limited to *Crypto* and *Giardia* with exception of BioSite Triage® which also detects *E. histolytica*

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

39

Rapid Cartridge Assays



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

40

2012-2013 Multi-site RCA Study


- WSLH participated in a CDC-sponsored study looking at the efficacy of the two major *Cryptosporidium* RCA's on the market at the time (Meridian Immunocard STAT! and Remel Xpect)
- Stools positive for *Cryptosporidium* by RCA at the clinical laboratory were submitted to WSLH; DFA (gold standard) and both RCA methods above were performed in-house

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

41

Evaluation of the Performance of Rapid Cartridge Assays for the Detection of *Cryptosporidium* from Stool

The TT, ALC3000, TA, Statens Serum Institut, University of Wisconsin-Labors, Madison, WI



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

42

2012-2013 Multi-site RCA Study

- 176 stool specimens evaluated at WSLH
- Meridian ImmunoCard STAT!
 - Sensitivity 94%
 - Specificity 73%
- Remel Xpect *Cryptosporidium*
 - Sensitivity 74%
 - Specificity 87%
- Similar data generated by the two other study sites

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

43

Single Target PCR Tests

- Many reference laboratories, state/regional public health laboratories and CDC have developed single target PCR assays to detect stool parasites
- Common targets are *Cryptosporidium* and *E. histolytica* (*histolytica*/*dispar*)
- High complexity, high costs to test and limitations of lab developed tests (LDT's) are some of the main reasons these tests did not become widely available

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

44

Multi-Target Stool Pathogen Tests

- Last two years has seen a major movement by numerous commercial test developers to get tests approved that will detect multiple stool pathogens simultaneously
- These rapid, multi-target tests have begun to re-shape the clinical testing and public health surveillance landscape
- Accurate and comprehensive testing is available without the need for traditional test methods on which current surveillance was defined/ built

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

45

Luminex xTAG GPP

- Both ASR and RUO kits available
- Besides detecting multiple viral and bacterial stool pathogens, has the ability to detect *Cryptosporidium*, *Giardia* and *E. histolytica*
- Bead-based technology; high throughput but labor intensive
- Next generation platform in development that will be more hands-off and efficient



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

46

BioFire FilmArray® GI Panel

- 22-target diagnostic GI panel; detects common bacterial, viral and the following parasitic agents:
 - *Cryptosporidium*
 - *Cyclospora*
 - *E. histolytica*
 - *Giardia*
- Multiplex PCR system
- Minimal hands-on time; results in approximately one hour; low throughput



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

47

BD Max™ Enteric Parasite Panel*

- Currently under development/ not launched*
- Will detect three common human stool parasites:
 - *Cryptosporidium*
 - *E. histolytica*
 - *Giardia*
- Flexible panel testing options; will complement enteric bacteria panel and enteric virus panel (also in development)



WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

48

Cryptosporidiosis Case Definition Changes - 2015

Updated CSTE case definition relies on laboratory diagnostic test used to diagnose infection.

Confirmed case= detection of **organisms or DNA** by:

- Direct fluorescent antibody (DFA).
- Polymerase Chain Reaction (PCR).
- Enzyme Immunoassay (EIA- **microplate only**).
- Light microscopy.

Probable case= **detection of antigen** by:

- Enzyme ImmunoAssay (EIA for antigen, microplate not specified).
- Immunochromatographic card test (i.e. ImmunoCard STAT!, some labs call these as EIAs).
- Rapid card test (some labs also call these EIAs).
- Unknown method.

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

49

Cryptosporidiosis Case Definition Changes- 2015

DPH survey of clinical laboratories to collect required information:

- Test kits used
- Testing protocols
- Referral practices

Development of classification algorithm to help local health departments with classification.

Benefits:

- More accurate confirmed vs. probable case numbers.
- Improved knowledge of how many labs are using rapid card tests.
 - Rapid card tests have a low positive predictive value and give many false positives.

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

50

Other Potential Impacts on Stool Parasite Surveillance

- Potential for improved surveillance due to detection of more parasites
 - Parasite testing may not be routinely ordered by clinicians in cases of GI illness
 - Expanded use of multi-target assays that include parasites will increase likelihood of their detection (*Cyclospora*, *Crypto*, *Giardia*); "syndromic" testing
- Significance of results in patients with multiple pathogens detected will have to be weighed

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

51

Molecular epidemiology- CryptoNet

- CDC launched molecular subtyping network for *Cryptosporidium* isolates.
- Differentiation/connection of outbreaks.
- Source tracking.
 - Linkage between cases and environmental source.
 - Identification of common zoonotic exposure.
- Improved epidemiologic understanding.
 - Geographic distribution
 - Common subtypes

***Subtyping can only be performed on specimens NOT fixed in formalin.**

CDC

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

52

Stool Specimen Submission- *Cryptosporidium*

- Stool specimens positive for *Cryptosporidium* are asked to be sent to WSLH for genotyping and surveillance testing
 - Fee exempt confirmation of positive specimens
 - Analysis of specimens unable to be confirmed or whose test results are in question
 - Genotyping performed for surveillance of clusters of illness
 - Can use Dunham Express courier service

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

53

Stool Specimen Submission- *E. histolytica/ dispar*

- If *E. histolytica/ dispar* is detected in a stool wet preparation or permanent stained smear, a clinician may inquire about the confirmation of the pathogenic *E. histolytica*
 - Some of the multi-target assays will differentiate
 - Reference laboratory performing the *E. histolytica* microplate EIA may differentiate
 - Submit unfixed stool (Raw, enteric culture transport, PVA) along with stool in formalin to WSLH for referral to CDC for PCR testing

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

54

Resources

Wisconsin State Laboratory of Hygiene

- Contact WSLH CDD Customer Service
- Dunham Express Courier**
- (800)236 – 7127
- Account 7271
- Next day delivery except on Sat/Sun
- Call WSLH customer service **(800)862-1013** during normal work hours or the WSLH pager service **(800)263-3280** after hours or weekends to discuss STAT testing

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

55

Resources

Wisconsin or national reference laboratories

- Many reference laboratories maintain the ability to perform comprehensive ova & parasite examinations of stool specimens
- Many may also be implementing multi-target PCR-based assays that will cover the common human parasitic pathogens that might be acquired and or detected locally

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

56

Resources

Centers for Disease Control and Prevention- Division of Parasitic Disease

- DPDx- Site maintained by the Division of Parasitic Diseases and Malaria (DPDM)
 - <http://www.cdc.gov/dpdx/>
- Can submit digital images to DPDx via email for telediagnosis; Fill out specimen submission form 50.34 (available from their web site or link from WSLH)
- Testing is fee exempt; STAT testing should be worked out directly with CDC contacts

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

57

Resources

- Information regarding which tests CDC offers can be acquired from the CDC web site:
 - <http://www.cdc.gov/laboratory/specimen-submission/list.html>
- If submitting directly to CDC, we ask that you fill out a WSLH req form (A) and fax it along with the CDC form 50.34 (if filled out) to WSLH (608-890-2548); All reports come back through WSLH; if already in the WSLH system, reporting will be expedited

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

58

CDC Test Directory

CDC Centers for Disease Control and Prevention
Submitting Specimens to CDC

Test Directory



CDC's Infectious Diseases Laboratories provides an online Test Directory that allows you to identify the right test for your needs. The searchable Test Directory includes an up-to-date list of clinical tests and procedures, information on specimen requirements, related information, test turnaround times, and other supplemental information. Access the directory here or while completing a Specimen Submission Form.

You may also download a copy of: www.cdc.gov/dpdx/pdfs/35811mainmenu-specimen-submission-cdc-10b.pdf

Effective December 31, 2014, an updated test directory is available. View the master list of changes: www.cdc.gov/dpdx/pdfs/35811mainmenu-specimen-submission-cdc-10b.pdf

Search
Narrow the results with a keyword, test title, test acronym, or point of contact:
[input field]

[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)
[Accutest](#)

Showing 32 of 32 items.

Test Name	Test Code
Accutest Molecular Detection	CDC-10471
Accutest Molecular Detection	CDC-10472
Accutest Molecular Detection	CDC-10473
Accutest Molecular Detection	CDC-10474
Accutest Molecular Detection	CDC-10475

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

59

Summary

- Traditional test methods are still utilized and effective for the detection of human stool parasites
- While not seen often, clinical laboratory parasitologists need to be able to recognize and identify human parasitic pathogens when they are present
- Numerous diagnostic resources are available via online (web-based), telediagnosis or via specimen referral.

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

60

Summary

- There are numerous commercial tests available for the detection of stool parasites in humans; **WSLH does not endorse any available test** but encourages laboratories to look at test performance, sensitivity, specificity and test validation when deciding to implement any such test
- Case definitions and surveillance are subject to change due to shifting landscape of diagnostic parasitology in WI and the U.S.

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

61

Contact Information

WSLH

Ann Valley- ann.valley@slh.wisc.edu

Tim Monson- timothy.monson@slh.wisc.edu

WSLH Customer Service- (800)862-1088

WDPH

Rachel Klos(608)267-7422

Communicable Epi Section(608)267-9003

CDC Division of Parasitic Diseases

(404) 718-4100/ dpd@cdc.gov

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

62

Thank You!

Thank
you

Google Images

Questions?



Google Images

WISCONSIN STATE LABORATORY OF HYGIENE - UNIVERSITY OF WISCONSIN

63